

AMINO ACIDS



The method for the analysis of amino acids by ion-exchange column and post-column derivatization was originally published in 1948. Its authors, Moore and Stein, were awarded the 1972 Nobel Prize in Chemistry.

Sixty years later, the analysis of a physiologic sample no longer takes two days and the resolution of the columns has been significantly improved. The method principles, however, remain the same.

One reason for this is the extraordinary insensitivity of the retention mechanism of the ion-exchange column against matrix effects, leading to a high reproducibility of retention times and quantification.

For amino acid analysis, Pickering offers the **PINNACLE PCX** and complete **Amino Acid Kits** containing eluents, reagents, diluents, columns and standards.

Two different types of kits are obtainable: **Sodium Kits** contain sodium citrate buffers and a sodium ion-exchange for the analysis of hydrolysates of peptides, proteins, collagens, foods and feeds.

Lithium Kits contain lithium citrate buffers and a lithium ion-exchange for the amino acid analysis of native samples. These are, for example, physiologic fluids (blood, urine) and cell tissues, but also food and drinks.

Derivatization is performed in one of two ways: either with TRIONE[®] (ninhydrin) or with o-phthalaldehyde (OPA) and Thiofluor[®] (2-mercaptoethanol derivative).

APPLICATION NOTE


Derivatization of Primary and Secondary Amino Acids with Ninhydrin

The most common reagent used for post-column derivatization of amino acids is Ninhydrin. Ninhydrin combined with primary amines generates a colored compound (Ruhemann's Purple; $\lambda_{\text{max}} = 570 \text{ nm}$; $\epsilon = 20000$). Hydrindantin serves as catalyst for the reaction, which is formed by reductive coupling from two molecules of ninhydrin.

A secondary amine reacts with Ninhydrin, even without the presence of hydrindantin, to form a yellow compound ($\lambda_{\text{max}} = 440 \text{ nm}$).


The detection limit (570 nm) for primary amino acids is approximately 30 pmol.

Is only a 1-channel-UV/VIS-detector available, both derivatives could also be detected with a slightly lower sensitivity at 500 nm.



Derivatization of Primary Amino Acids with *o*-Phthalaldehyde/Thiofluor®

Primary amines can be converted with *o*-phthalaldehyde (OPA) and Thiofluor (2-mercaptoethanol derivative) in an alkaline environment (pH 9-10) to a fluorescent isoindole derivative. The detection limit of with OPA derivatized amino acids is about 5 pmol.



HPLC Conditions and Derivatization Parameters

Derivatization with Ninhydrin (TRIONE®)

| HPLC | |
|-----------------------------------|--|
| Operational Mode | quarternary gradient |
| Eluant | Lithium or sodium buffer |
| Degassing | Helium or vacuum degassed |
| HPLC Column | Lithium or sodium ion-exchange columns with GARD |
| Flow Rate | 0.3 – 0.6 ml/min (depends on the method) |
| Injection Volume | Up to 100 µL |
| Post-Column Derivatization | |
| Pinnacle PCX | Single-pump |
| Column Oven | 34 – 70 °C (gradient depends on the method) |
| Reactor Volume | 500 µL |
| Reactor Temperature | 130 °C |
| Reagent | TRIONE® (Ninhydrin/Hydrindantin) |
| Reagent Flow | 0.3 – 0.55 mL/min (depends on the method) |
| Detection | |
| Measurement | UV/VIS-Detector |
| UV/VIS | 570 nm; primary amino acids |
| UV/VIS | 440 nm; secondary amino acids |
| Flowcell | Analytical; pressure stable up to 7 bar |

APPLICATION NOTE

Derivatization with o-Phthalaldehyde/Thiofluor®

| HPLC | |
|-----------------------------------|---|
| Operational mode | quarternary gradient |
| Eluent | Lithium or sodium buffer |
| Degassing | Helium or vacuum degassed |
| HPLC-column | Lithium or sodium ion-exchange column with GARD |
| Flow Rate | Depends on the method (0.3 -0.6 ml/min) |
| Injection Volume | Up to 100 µL |
| Post-Column Derivatization | |
| Pinnacle PCX | Single-pump |
| Column Oven | 34 – 70 °C (gradient depends on the method) |
| Reactor Volume | 150 µL |
| Reactor Temperature | 45 °C |
| Reagent | o-Phthalaldehyde, Thiofluor® |
| Reagent Flow | 0.3 mL/min |
| Detection | |
| Detection Type | Fluorescence detection |
| Excitation Wavelength | 330 nm |
| Emission Wavelength | 465 nm |
| Flowcell | Analytical; pressure stable up to 7 bar |

Caution: Extreme pH-range!

As a consequence of the alkaline regenerant (pH 13), components made from Vespel may not be present in the HPLC system. These must be exchanged for components made from pH-inert materials (Tefzel or PEEK). For advice, please contact your LC-representative or consult your technical manual.

An inert version (Titanium, PEEK) of the HPLC system is not necessary, a piston seal wash for the pump heads, however, is recommended.

To avoid corrosion of the system and contamination of the ion-exchange column with metal ions, passivating is recommended; this is especially valid for older systems. For advice, please contact your LC-representative or consult your technical manual.

Autosampler

To achieve reproducible retention times in amino acid analysis with ion-exchange columns, it is important to keep to exact time intervals between injections. Should only a manual injector be available, then the injection times should be controlled with a laboratory timer.

APPLICATION NOTE

Amino Acid Analysis Using Sodium Ion-exchange Columns

For amino acid analysis of hydrolysates of peptides, proteins, collagens, foods and feeds we offer **Sodium Kits** with sodium citrate buffers and an ion-exchange column (type "Standard" or "High Efficiency") in sodium form.

The following kits are available:

For Pinnacle PCX Instrument (with the possibility to use a temperature gradient):

| | | |
|---|---------------|------------------|
| Protein-Hydrolysates (33-Minutes Chromatogram) | Trione T200 * | 0352-0058 |
| | Trione T100C* | 0352-0057 |
| | OPA* | 0352-0059 |
| Collagen hydrolysates (33-Minutes Chromatogram) | Trione T200 | 0352-0062 |
| | Trione T100C | 0352-0061 |
| | OPA | 0352-0063 |
| Oxidized feed hydrolysates (33-Minutes Chromatogram) | TRIONE T200 | 0352-0021 |
| | TRIONE T100C | 0352-0020 |
| | OPA | 0352-0022 |

* Trione T200: shelf life: 12 months; Trione T100C: shelf life: 3 months; OPA: o-Phthalaldehyde

APPLICATION NOTE

Method descriptions

Please use the HPLC conditions and derivatization parameters on page 3ff if there are no other parameters given.

Method: Protein- and Collagen-Hydrolysates

Protein- and Collagen-Hydrolysates: 33-minutes-chromatogram:

"High Efficiency" Ion-exchange column (1154110T); 4.6 x 110 mm; Na⁺ Form

Pickering-Kit 0352-0058, 0352-0057 or 0352-0059 for protein hydrolysates

Pickering-Kit 0352-0062, 0352-0061 or 0352-0063 for collagen hydrolysates

Parameters:

| | | |
|-----------------------------|----------|---|
| Flow rates: | Eluant: | 0.6ml/min |
| | Reagent: | 0.3ml/min |
| Column initial temperature: | | 46°C for protein hydrolysates 42°C for collagen hydrolysates |

HPLC-Program:

| Step | Time | Interval | Na315 % | Na425 % | Na640 % | RG011 % | Comment |
|--------------|------|----------|---------|---------|---------|---------|------------------|
| <i>Equil</i> | | | 100 | | | 0 | Equilibration |
| 1 | 0 | 0 | 100 | 0 | 0 | 0 | Inject 10 µL |
| 2 | 4.0 | 4 | 100 | 0 | 0 | 0 | Isocratic |
| 3 | 15.0 | 11 | 0 | 100 | 0 | 0 | Linear gradient |
| 4 | 16.0 | 1.0 | 0 | 0 | 100 | 0 | Linear gradient |
| 5 | 31.0 | 15.0 | 0 | 0 | 100 | 0 | Isocratic |
| 6 | 31.1 | 0.1 | 0 | 0 | 0 | 100 | Step change |
| 7 | 33.0 | 1.9 | 0 | 0 | 0 | 100 | Clean-out |
| 8 | 33.1 | 0.1 | 100 | 0 | 0 | 0 | Re-equilibration |
| 9 | 40 | 6.9 | 100 | 0 | 0 | 0 | Re-equilibration |

Run time 33 min

Equilibration time 7 min

Column oven program for protein hydrolysates:


| Time | Temperature, °C |
|------|-----------------|
| 0 | 46 |
| 4 | 46 |
| 9 | 70 |
| 32 | 70 |
| 33 | 46 |

Column oven program for collagen hydrolysates:


| Time | Temperature, °C |
|------|-----------------|
| 0 | 42 |
| 4 | 42 |
| 9 | 70 |
| 32 | 70 |
| 33 | 42 |

APPLICATION NOTE

Chromatogram of an Amino Acid Standard for Protein Hydrolysates
 "High Efficiency" Ion-exchange Column, 4.6 x 110 mm, Na⁺ Form
 Sample: Protein Hydrolysate Standard, 0.25 µmole/mL (Catalog No.: 012506H)



Chromatogram of an Amino Acid Standard for Collagen Hydrolysates
 "High Efficiency" Ion-exchange Column, 4.6 mm x 110 mm, Na⁺ Form
 Sample: Collagen Hydrolysate Standard, 0.25 µmole/mL (Catalog No.: 012506C)



APPLICATION NOTE

Method: Oxidized Feed Hydrolysates

Oxidized Feed Hydrolysates: 33-minutes-chromatogram

High Efficiency" Ion-Exchange-column (1154110T); 4.6 x 110 mm, Na⁺-Form

Pickering-Kit 0352-0021, 0352-0020 or 0352-0022 for oxidized feed samples

Parameters:

Flow rates: Eluant: 0.6 ml/min
 Reagent: 0.5 ml/min
 Column initial temperature: 55°C

HPLC-Program:

| Step | Time | Interval | Na270 % | Na425 % | Na640 % | RG011 % | Comment |
|--------|------|----------|---------|---------|---------|---------|------------------|
| Equil. | | | 100 | 0 | 0 | 0 | Equilibration |
| 1 | 0 | 0 | 100 | 0 | 0 | 0 | Inject 10 µL |
| 2 | 4.0 | 4 | 100 | 0 | 0 | 0 | Isocratic |
| 3 | 15.0 | 11 | 0 | 100 | 0 | 0 | Linear gradient |
| 4 | 16.0 | 1 | 0 | 0 | 100 | 0 | Linear gradient |
| 5 | 31.0 | 15 | 0 | 0 | 100 | 0 | Isocratic |
| 6 | 31.1 | 0.1 | 0 | 0 | 0 | 100 | Step change |
| 7 | 33.0 | 1.9 | 0 | 0 | 0 | 100 | Cleanout |
| 8 | 33.1 | 0.1 | 100 | 0 | 0 | 0 | Re-equilibration |
| 9 | 40 | 6.9 | 100 | 0 | 0 | 0 | Re-equilibration |

Run time 33 min


Equilibration time 7 min

Column Oven Program:

| Step | Time | Temperature, °C |
|------|------|-----------------|
| 1 | 0 | 55 |
| 2 | 12 | 55 |
| 3 | 17 | 70 |
| 4 | 32 | 70 |
| 5 | 33 | 55 |

APPLICATION NOTE

Chromatogram of an Amino Acid Standard for oxidized feed hydrolysates
 "High Efficiency" Ion-Exchange Column, 4.6 mm x 110 mm, Na⁺-Form (1154110T)
 Sample: Oxidized Feed Hydrolysate Standard, 0.25 µmole/ML (Catalog No. 1700-0155)



Several aminoacids of interest in feed can be analyzed in a significant shorter runtime:


➤ **Method for accelerated Lysine analysis (9min)**

Parameters:

Flow rates: Eluant: 0.6ml/min
 Reagent: 0.5ml/min
 Column temperature: 60°C (no temperature gradient necessary)

HPLC-Program:

| Time , min | % Na640 | % RG011 |
|------------|---------|---------|
| 0 | 100 | 0 |
| 6 | 100 | 0 |
| 6.1 | 0 | 100 |
| 9 | 0 | 100 |
| 9.1 | 100 | 0 |
| 15 | 100 | 0 |



Run time: 9 min
 Equilibration time: 6 min

APPLICATION NOTE

➤ **Method for accelerated analysis of Cysteic Acid, Methionine Sulfone and Lysine (16 min)**


Parameters:

Flow rates: Eluant: 0.6ml/min
 Reagent: 0.5ml/min
 Column temperature: 55°C (no temperature gradient necessary)

HPLC-Program:

| Time , min | % Na270 | % Na640 | % RG011 |
|------------|---------|---------|---------|
| 0 | 100 | 0 | 0 |
| 5 | 100 | 0 | 0 |
| 5.1 | 0 | 100 | 0 |
| 12 | 0 | 100 | 0 |
| 12.1 | 0 | 0 | 100 |
| 16 | 0 | 0 | 100 |
| 16.1 | 100 | 0 | 0 |
| 23 | 100 | 0 | 0 |

Run time 16 min
 Equilibration time 7 min



APPLICATION NOTE

Amino Acid Analysis Using Lithium Ion-Exchange Columns

For amino acid analysis of native samples, e.g. blood, urine and tissues, and also foods and beverages, we offer **Lithium Kits** with lithium citrate buffers and an ion-exchange column in lithium form.

The following kit is available:

For Pinnacle PCX with the possibility to use a temperature gradient:

| Method | | Order Number |
|--|---------------|--------------|
| Native Samples (70-minutes-chromatogram)* | TRIONE® T200 | 0352-0007 |
| | TRIONE® T100C | 0352-0006 |
| | OPA | 0352-0008 |

* It is not possible to use Norleucine as internal standard; we suggest to use for example Glucosaminic Acid

Method: Physiologic Fluids and Native Samples

Physiological Fluids: 70-minutes-chromatogram

"High Efficiency" Ion-exchange column (0354675T); 4.6 x 75 mm; Li⁺ Form

Pickering-Kit 0352-0006, 0352-0007 or 0352-0008 for physiological fluids

Parameters:

Flow Rate: Eluant: 0.55 ml/min
Column initial temperature: 34°C

HPLC program:

| Step | Time | 1700-1125 % | Li365 % | Li375 % | RG003 % | Comment |
|------|------|-------------|---------|---------|---------|------------------|
| 0 | 0 | 100 | 0 | 0 | 0 | Inject |
| 1 | 10 | 100 | 0 | 0 | 0 | Isocratic |
| 2 | 19 | 40 | 60 | 0 | 0 | Linear Gradient |
| 3 | 32 | 0 | 100 | 0 | 0 | Linear Gradient |
| 4 | 43 | 0 | 100 | 0 | 0 | Isocratic |
| 5 | 43.1 | 0 | 0 | 100 | 0 | Step Change |
| 6 | 57 | 0 | 0 | 100 | 0 | Isocratic |
| 7 | 57.1 | 0 | 0 | 70 | 30 | Step Change |
| 8 | 72 | 0 | 0 | 70 | 30 | Isocratic |
| 9 | 72.1 | 100 | 0 | 0 | 0 | Step Change |
| 10 | 88 | 100 | 0 | 0 | 0 | Re-equilibration |

Run time: 72 min


Equilibration time: 16 min

APPLICATION NOTE

Column oven program:

| Time, min | Temperature, °C |
|-----------|-----------------|
| 0 | 34 |
| 6 | 34 |
| 17 | 65 |
| 25 | 70 |
| 70 | 70 |
| 71 | 34 |

Chromatogram of an Amino Acid Standard for Native Samples
 "High Efficiency" Ion-exchange column, 4.6 x 75 mm, Li⁺ Form
 Sample: Lithium Calibration Standard, 0.25 µmole/mL (Catalog No. 1700-0170)



APPLICATION NOTE

Chemicals and Columns

Post-Column Derivatization Unit

| Order Number | Description |
|--------------|---|
| 1153-1022 | PINNACLE PCX; single-pump, 500 µL reactor for TRIONE® |
| 1153-1012 | PINNACLE PCX; single-pump, 150 µL reactor for OPA |

Kits for the Analysis of Hydrolysates

For all mentioned kits below, it is necessary to use a temperature gradient of the column:

| Order Number | Description | |
|---|---|--------------|
| 0352-0057 | 30-Min High Efficiency Kit for Protein Hydrolysates | |
| 1154110T | Na ⁺ Ion-Exchange Column | 4.6 x 110 mm |
| 1700-3102 | Kit, GARD Holder & Cation GARD | 2 units |
| Na315 | Na ⁺ Eluent, pH 3.15 | 4 x 950 mL |
| Na425 | Na ⁺ Eluent, pH 4.25 | 4 x 950 mL |
| Na640 | Na ⁺ Eluent, pH 6.40 | 4 x 950 mL |
| RG011 | Na ⁺ Column Regenerant | 950 mL |
| T100C | TRIONE® Ninhydrin Reagent, ready to use (max 3-month shelf life from date of production) | 4 x 950 mL |
| Na220 | Na ⁺ Diluent, pH 2.20 | 4 x 250 mL |
| 012506H | Na ⁺ Calibration Standard for protein hydrolysate, 0.25 µmol/mL | 5 mL |
| | | |
| 0352-0061 | 30-Min High Efficiency Kit for Collagen Hydrolysates | |
| Kit identical to 0352-0057, but replaces the standard with: | | |
| 012506C | Na ⁺ Calibration Standard for collagen hydrolysate, 0.25 µmol/mL | 5 mL |
| | | |
| 0352-0058 | Kit identical to 0352-0057 with T200 replacing T100C | |
| 0352-0062 | Kit identical to 0352-0061 with T200 replacing T100C | |
| T200 | TRIONE® Two-part Ninhydrin Reagent, to prepare | 4x 900 mL |
| | | |
| 0352-0059 | Kit identical to 0352-0057 with items below replacing T100C | |
| 0352-0063 | Kit identical to 0352-0061 with items below replacing T100C | |
| O120 | o-Phthalaldehyde (OPA), | 5 g |
| OD104 | OPA Diluent, | 4 x 950 mL |
| 3700-2000 | Thiofluor™ (2 each per kit) | 10 g |

APPLICATION NOTE


| | | |
|------------------|---|---------------|
| 0352-0021 | 33-Min High Efficiency Kit for oxidized feed hydrolysates | |
| 1154110T | Na ⁺ Ion-Exchange-column | 4.6 x 110 mm |
| 1700-3102 | Kit, Gard Holder & Cation GARD | 2 pieces/pack |
| Na270 | Na ⁺ Eluant, pH 2.70 | 4 x 950 mL |
| Na425 | Na ⁺ Eluant, pH 4.25 | 4 x 950 mL |
| Na640 | Na ⁺ Eluant, pH 6.40 | 4 x 950 mL |
| RG011 | Na ⁺ Column Regenerant | 950 mL |
| T100C | TRIONE® Ninhydrin Reagent, ready to use (max 3-month shelf life from date of production) | 4 x 950 mL |
| Na220 | Na ⁺ Diluent, pH 2.20 | 4 x 250 mL |
| 1700-0155 | Na ⁺ calibration standard for oxidized feed hydrolysates; 0.25 µmol/mL | 5 mL |
| | | |
| 0352-0020 | Kit identical to 0352-0021 with T100C replacing T200 | |
| T200 | TRIONE® Two-part Ninhydrin Reagent, to prepare | 4x 900 mL |
| | | |
| 0352-0022 | Kit identical to 0352-0021 with items below replacing T100C | |
| O120 | o-Phthalaldehyde (OPA), | 5 g |
| OD104 | OPA Diluent, | 4 x 950 mL |
| 3700-2000 | Thiofluor™ (2 each per kit) | 10 g |

Kits for the Analysis of Physiologic Fluids and Native Samples

| Order Number | Description | |
|------------------|---|-------------|
| 0352-0006 | 70-min Physiologic Fluid/Native Sample Analysis Kit | |
| 0354675T | Li+ Ion-Exchange Column & 1700-0070 AA test mix | 4.6 x 75 mm |
| 1700-3102 | Kit, GARD Holder & Cation GARD | 2 units |
| 1700-1125 | Li+ Eluant, pH 2.80 | 4 x 950 mL |
| Li365 | Li+ Eluant, pH 3.65 | 4 x 950 mL |
| Li375 | Li+ Eluant, pH 3.75 | 4 x 950 mL |
| RG003 | Li+ Column Regenerant | 950 mL |
| T100C | TRIONE® Ninhydrin Reagent, ready to use | 4 x 950 mL |
| Li220 | Li+ Diluent, pH 2.20 | 4 x 250 mL |
| SP100 | SERAPREP™ | 250 mL |
| UP100 | URIPREP™ | 250 mL |
| 1700-0170 | Li+ Calibration Standard, without Norleucin and alpha-Amino-β-guanidinopropionic acid, Physiologic fluid, 0.25 μmole/mL | 5ml |
| | | |
| 0352-0007 | Kit identical to 0352-0006 with T200 replacing T100C | |
| T200 | TRIONE® Two-part Ninhydrin Reagent, to prepare | 4 x 900 mL |
| | | |
| 0352-0008 | Kit identical to 0352-0006 with items below replacing T100C | |
| O120 | o-Phthalaldehyde (OPA) | 5 g |
| OD104 | OPA Diluent | 4 x 950 mL |
| 3700-2000 | Thiofluor™ (2 each per kit) | 10 g |

Caution-Exchange Columns for Amino Acid Analysis

Each amino acid column is tested by applying an amino acid calibration standard and running it according to a specific gradient elution protocol. After a column passes this rigorous QC test, it receives a serial number and is packaged with its test chromatogram. Pickering Laboratories specializes in the manufacture of cation-exchange columns for amino acid analysis. These are described in more detail below. Although pre-column derivatization followed by reversed-phase chromatography is popular with some HPLC users, its successful application is limited to protein hydrolysates due to matrix interference with native samples.



Without extensive pre-treatment, the carbohydrates and organic acids in wine, the fats in tobacco, and the macromolecules in human serum will cause the amino acids to elute at a different time for each unique matrix.

In the reversed-phase environment all species in solution compete simultaneously for residency in the stationary phase. The retentiveness of any particular species will be influenced by the concentration of all other solutes present.

In contrast, ion-exchange chromatography followed by post-column derivatization is intrinsically more rugged and repeatable than pre-column analysis. As chromatography and derivatization represent two separate and distinct events, each can be optimized independently.

The retention mechanism in ion-exchange chromatography is almost completely matrix-insensitive. Since the sulfonated divinylbenzene polymer comprising the stationary phase has a high ion-exchange capacity, the positively charged amino acids exhibit a strong affinity for these fixed negative sites. Consequently, all amino acids, except for the very acidic, are fixed in a narrow band in the guard column, while the sample matrix moves on. Thus the same program used by one laboratory for the analysis of wine can be employed by another to analyze human urine.

The most important difference between ion-exchange resins and reversed-phase silicas is in chemical selectivity. Bonded reversed-phase silicas are manufactured to exhibit monotypical chromatographic behavior, e.g., partitioning.

In contrast ion-exchange resins are polytypical. Separation is based not only upon ion-exchange, but also partitioning, adsorption, ion exclusion, and more. Because of the involvement of so many retention mechanisms, a single change in any operational parameter—cation concentration, pH, flow rate, column temperature—can result in multiple changes in peak position. This flexibility is particularly advantageous, for example, when developing a method to optimize the separation of a few amino acids of interest at the expense of the others, as in the PKU method and other rapid screens.

Pickering also offers also a new guard column, the cation exchange column **GARD™**:

The new GARD™ guard column is recommended to protect the HPLC column from contamination by the matrix. This will avoid significant costs. The GARD™ prolongs the life time of your HPLC column without band spreading or additional pressure. The GARD™ precolumn is very easy to use and also more cost-effective for your laboratory and can be used for almost all cation-exchange applications.

Cation-Exchange Columns for Analysis of Hydrolysates

- Constant Column Temperature or Temperature Gradient
- "High Efficiency" 30-Minute (with temperature gradient) or 55-Minute (without temperature gradient) Analysis

| Catalog No. | Description |
|-------------|--|
| 1154110T | "High efficiency" Analytical Column (30 minutes) for the analysis of protein and collagen hydrolysates or sulfur containing amino acids in oxidized feed hydrolysates; cation-exchange, Na ⁺ form, 4.6 x 110 mm (no guard column available) |
| 1193250 | Standard Analytical Column (60 minutes) for the analysis of protein and collagen hydrolysates; cation-exchange, Na ⁺ form, 3 x 250 mm |
| 1154150T | "High efficiency" Analytical Column (55 minutes) for the analysis of protein and collagen hydrolysates; cation-exchange, Na ⁺ form, 4 x 150 mm |
| 1700-3102 | GARD: consists of holder and 2 columns |

Cation-Exchange Columns for the Analysis of Native Resp. Physiologic Samples

- Constant Column Temperature or Temperature Gradient
- "High Efficiency" 70-Minutes Analysis or "High Efficiency" 95-Minute Analysis (both with temperature gradient)

| Order Number | Description |
|--------------|---|
| 0354675T | "High Efficiency" Analytical Column (70 minutes) for the analysis of native resp. physiologic samples; cation-exchange, Li ⁺ form, 4.6 x 75 mm (no guard column available) |
| 0354100T | "High Efficiency" Analytical Column (95 minutes or 120 minutes) for the analysis of native resp. physiologic samples; cation-exchange, Li ⁺ form, 4.0 x 100 mm |
| 1700-3102 | GARD: Kit consists of holer and 2 columns |

TRIONE® Ninhydrin Reagent



- 1) Stored at room temperature
- 2) Superb signal/noise-ratio
- 3) Superb signal/noise-ratio, thereby high detection sensitivity

TRIONE is a Ninhydrin reagent especially formulated for amino acid analysis. It is so stable that it does not require refrigeration, either in shipment, storage, or in the reservoir. Quantitation is consistent from the first to the last milliliter, eliminating waste. The high signal-to-noise ratio of TRIONE, when compared to DMSO-containing reagents, permits detection sensitivity to be increased with minimum increase in background noise - a feature particularly appreciated at sample amounts of <50 picomoles.

TRIONE Ninhydrin Reagent is a patented formulation containing Ninhydrin, Hydrindantin (reduced Ninhydrin), a Lithium Acetate buffer, and sulfolane, a water-miscible organic solvent. The solvent is necessary to maintain the solubility of both the Hydrindantin and the primary amine product, Ruhemann's Purple. The buffer is required because the reaction is pH dependent. The active ingredients - Ninhydrin and Hydrindantin - are required for proper

development of secondary and primary amines, respectively.

Two preparations are available to suit your usage and storage requirements:

T100C

- 1) Pour into your reservoir and use; the ultimate in convenience with a minimum of handling
- 2) Three-month shelf life


T200

- 1) Combine two solutions, swirl and use
- 2) Twelve-month shelf life before mixing; three months in the reservoir

| Order Number | Description |
|--------------|--|
| T100C | TRIONE Ninhydrin Reagent, 4 x 950 mL |
| T200 | TRIONE Ninhydrin, two solutions for mixing, gives 4 x 950 mL |

APPLICATION NOTE

o-Phthalaldehyde Reagent



Primary amines such as amino acids form highly fluorescent compounds when reacted with o-phthalaldehyde (OPA) and a mercaptan under basic conditions. At a pH >9 and ambient temperature, reaction is generally complete within 1 - 30 seconds. The products of this reaction, 1-alkyl-2-alkylthio-substituted isoindoles, exhibit optimal excitation at 330 nm and maximal emission at 465 nm.

For an oxygen-sensitive reagent like OPA to remain stable for days instead of hours, it is important to start with the purest and most stable ingredients available, and to store and use the reagent under anaerobic conditions.

Using the chemicals described below, a long-lasting (up to ten days) OPA reagent for post-column derivatization of primary amines can be prepared. Each chemical is accompanied by clear instructions for formulating the reagent in your laboratory within minutes.

For the preparation of o-phthalaldehyde reagent PICKERING offers:

- 1) Ultra-pure, crystalline OPA
- 2) Borate diluent; free of heavy metals, particulates and amines
- 3) Thiofluor; crystalline substitute for 2-mercaptoethanol

o-Phthalaldehyde

| Order Number | Description |
|--------------|---|
| O120 | o-Phthalaldehyde (OPA), "Chromatographic Grade™", 5 g |

o-Phthalaldehyde Diluent

For the amino acid analysis Pickering offers a sodium borate buffer with pH 10.4 capable to buffer acidic eluents (pH 2).

| Order Number | Description |
|--------------|---|
| OD104 | OPA Diluent, "Chromatographic Grade™", 4 x 950 mL |

Thiofluor[®]

Pickering's Thiofluor, a solid, nearly odorless nucleophile, is a superior substitute for 2-mercaptoethanol in the preparation of OPA reagents. It forms a more stable and longer-lasting fluorophore with OPA than does 2-mercaptoethanol, while possessing the same fluorescence properties.

Unlike the volatile 2-mercaptoethanol, Thiofluor will not migrate through the gas manifold and regulator of the OPA reagent pressurization system.

| Order Number | Description |
|--------------|--|
| 3700-2000 | Thiofluor [®] , "Chromatographic Grade [™] ", 10 g |

Eluent Buffers and Regenerants

- 1) "Chromatographic Grade[™]"
- 2) Guaranteed pH stability at ambient conditions
- 3) Negligible amine contamination for optimum signal-to-noise characteristics
- 4) Cost-effective: long shelf life; can be used to the last milliliter
- 5) Consistent elution profiles, bottle-to-bottle, lot-to-lot
- 6) No refrigeration required

There are two dominant mechanisms for effecting separation in ion-exchange chromatography. One mechanism is *titration*, or increasing the pH. As pH increases, the amino acids go from a positive charge to a neutral state, and so are released from the fixed anionic sites in the resin (stationary phase).

The other mechanism is *competitive ion-exchange*. The positively charged amino acids show a higher affinity for the ion-exchange sites than do the eluting cations (H^+/Li^+ or Na^+) and so they are bound strongly to the anionic sites. They are eluted by the continuous flow of cation or by increasing cation concentration developed by gradient formation. The earliest protocols for ion-exchange of amino acids were based on step gradients—*isocratic* procedures that step through a set of 3–6 solutions of varying pH and cation concentration. Both titration and elution adjustments are made with each eluant step. This technique persists today in the design of dedicated amino acid analyzers.

The modern liquid chromatographic (HPLC) techniques for ion-exchange use three solutions and continuous gradients. Gradient capability allows for easier methods development and results in flatter baselines—especially important at high sensitivity. As with dedicated analyzers, the eluants are either sodium- or lithium-based. The need for two cation systems derives from the two fundamental types of samples encountered in amino acid analysis: hydrolyzed and non-hydrolyzed ("native") samples. Hydrolyzing the sample results in a smaller array of amines (typically up to 21) than is present in native samples (typically 40–60). The wider range of amines in native samples requires a more discriminating or weaker eluant to fully resolve them. A lithium ion-based eluant provides the necessary resolution.

If the sample is hydrolyzed, a sodium system can be used for a shorter run time. Both Na^+ and Li^+ eluants for gradient HPLC are designed in the same way. The first, low-pH eluant is a buffer with low cation concentration. The second eluant has no buffering capacity and a high cation concentration. The third, high-pH eluant (usually referred to as the column regenerant) is based on the cation hydroxide and is of comparable cation

APPLICATION NOTE

concentration to the first eluant. These solutions allow for titration relatively independent of changes in eluant strength, or allow for increasing eluant strength independent of pH.

Pickering's Lithium and Sodium eluants are not sensitive to oxidation and do not need refrigeration, either in storage or use. Degassing is not required. However, they should be protected from air to prevent contamination. Ambient air actually contains amines and amino acids that will dissolve in the low-pH eluants and will appear in the chromatograms. All buffers are packaged in cases of four 950 mL polyethylene bottles. The cation hydroxide regenerants are packaged in single 950 mL bottles because of the small volumes consumed during each analysis. Sample diluents are packaged in cases of four 250 mL borosilicate glass bottles.

| Order Number | Description |
|--|--|
| Sodium buffers for hydrolysates | |
| Na270 | Sodium Eluent, pH 2.70, 4 x 950 mL |
| Na315 | Sodium Eluent, pH 3.15, 4 x 950 mL |
| Na328 | Sodium Eluent, pH 3.28, 4 x 950 mL |
| Na425 | Sodium Eluent, pH 4.25, 4 x 950 mL |
| Na640 | Sodium Eluent, pH 6.40, 4 x 950 mL |
| Na740 | Sodium Eluent, pH 7.40, 4 x 950 mL |
| 1700-0112* | Sodium Eluent, pH 3.15, 4 x 950 mL, with 5 % sulfolane |
| 1700-0114* | Sodium Eluent, pH 3.15, 4 x 950 mL, with 2.5 % sulfolane |
| RG011 | Column regenerant, Na ⁺ Form, 950 mL |
| * for use with 1154150 with Serial No. above 1314 | |
| Lithium buffer for native and physiologic samples | |
| Li275 | Lithium Eluent, pH 2.75, 4 x 950 mL |
| 1700-1125 | Lithium Eluent, pH 2.80, 4 x 950 mL |
| Li280 | Lithium Eluent, pH 2.80, 4 x 950 mL |
| Li292 | Lithium Eluent, pH 2.92, 4 x 950 mL |
| Li357 | Lithium Eluent, pH 3.57, 4 x 950 mL |
| Li365 | Lithium Eluent, pH 3.65, 4 x 950 mL |
| Li375 | Lithium Eluent, pH 3.75, 4 x 950 mL |
| Li750 | Lithium Eluent, pH 7.50, 4 x 950 mL |
| RG003 | Column regenerant, Li ⁺ Form, 950 mL |

Sample Diluents

Use of these sample diluents is essential to ensuring reproducibility from injection to injection. They establish a uniform pH and ion concentration at the outset, regardless of the source and pre-treatment of the sample. The sample is maintained buffered and at optimum pH for sample storage and analysis. Use Na220 for hydrolysates and Li220 for native amino acid samples.

| Order Number | Description |
|--------------|--|
| Na220 | Sodium Diluent, pH 2.20, 4 x 250 mL, "Chromatographic Grade™" |
| Li220 | Lithium Diluent, pH 2.36, 4 x 250 mL, "Chromatographic Grade™" |

SERAPREP® und URIPREP®

- 1) "Chromatographic grade™"
- 2) For the preparation of native samples
- 3) Ensures reproducibility of early eluting peaks
- 4) Accommodates sample buffering capacity
- 5) Little to no pH adjustments
- 6) Mix, spin, filter, inject



Native amino acids are those found "free" in samples such as serum, urine and other physiologic fluids, plant extracts, foods and beverages. Although preparation of these samples for amino acid analysis is much simpler and less time-consuming than protein hydrolysis, control of pH and normality, and removal of soluble protein are critical factors which can affect the chromatography.

The early-eluting amino acids - taurine, urea, aspartic acid, threonine, serine, etc. - are particularly sensitive to pH and normality. Accordingly, the samples must be held to a narrow pH range between 2.1 and 2.5, and at the proper Lithium ion concentration to ensure reproducibility in the early part of the chromatogram. The later-eluting compounds are more tolerant of initial sample conditions, and their retention times are not as likely to be affected.

SERAPREP® and URIPREP® replace commonly used protein precipitation reagents such as acetonitrile, trichloroacetic acid and picric acid, and eliminate the need for dialysis, ultrafiltration, and repeated centrifugation steps, followed by pH adjustment.

Use SERAPREP® for preparing serum and URIPREP® for preparing urine and other samples with low buffering capacity, such as orange juice or wine. The efficiency of protein precipitation and need for post-centrifugation pH adjustment of the sample determine which reagent is best for your particular sample.

APPLICATION NOTE

In a microcentrifuge tube, thoroughly mix equal portions of sample and SERAPREP® or URIPREP®.

1. Let stand for 5 minutes. Centrifuge the mixture at 13.000 rpm for 5 minutes. Check the supernatant pH to ensure that the range is pH 2.3 ±0.2. Adjust the initial mixing ratio as necessary.
2. Filter the supernatant with a syringe filter (0.2 or 0.45 µm). The filtrate is ready to be injected into an auto-sampler vial for amino acid analysis.
3. If further dilution is needed, use Li 220 to adjust the concentration of analyte.

| Order Number | Description |
|--------------|---|
| SP100 | SERAPREP®, Sample Preparation for Serum, 250 mL |
| UP100 | URIPREP®, Sample Preparation for Urine, 250 mL |

Amino Acid test Mixture (1700-0070)

L-Cysteic Acid
L-Threonine
L-Serine

This test mixture is in 0.01 N HCl at a concentration of 0.25 µmol/mL. The three amino acids are eluted isocratically with the first eluant (Na328, Na315, Li275 or Li280). Since cysteic acid is not retained by the stationary phase, its retention time may be used to calculate the system void volume. The resolution between threonine and serine indicates the degree of band-spreading due to system leakages or excessive dead volumes.

The Test Mixture is a convenient tool for rapid, routine monitoring of the performance of either a sodium or lithium amino acid analysis system.

Test mixture should not be used for quantitation

Amino Acid Standards

- 1) Quantitative
- 2) Chromatographically-pure starting components
- 3) Each lot tested by amino acid analysis

Pickering's amino acid standards have a reputation worldwide for quality and reliability in all post-column systems and methods. Six amino acid mixtures are available for a variety of applications. Except for the three-component test mixture (1700-0700) the standards are in 5 mL vials, in an appropriate citrate buffer.

Although they are stored frozen at the factory, the Pickering calibration standards will remain stable when shipped at ambient temperatures. Upon receipt, however, it is important to place them into a freezer immediately and store frozen until ready for use.

APPLICATION NOTE



| CALIBRATION STANDARDS FOR AMINO ACID ANALYSIS | | | | | | | | | |
|---|------------|------------|---------|---------|-----------|----------|---------|-----------|-----------|
| CONSTITUENTS | 1700-0180* | 1700-0175* | 011006P | 012006P | 1700-0150 | 012506C | 012506H | 1700-0155 | 1700-0170 |
| Beta-Alanine | * | | * | * | | | | | * |
| Alanine | * | | * | * | | * | * | * | * |
| D,L-a-Amino-adipic acid | * | | * | * | | | | | * |
| Gamma-Amino butyric acid | | * | * | * | | | | | * |
| Alpha-Amino-n-butyric acid | * | | * | * | | | | | * |
| D,L,b-Amino-i-butyric acid | * | | * | * | | | | | * |
| Alpha-Amino-Beta-guanidinopropionic acid | | | * | * | | | | | |
| Ammonia | | * | * | * | | * | * | * | * |
| Anserine | | * | * | * | | | | | * |
| Arginine | | * | * | * | | * | * | * | * |
| Asparagine | * | | * | * | | | | | * |
| Aspartic acid | * | | * | * | | * | * | * | * |
| Carnosine | | * | * | * | | | | | * |
| Citrulline | * | | * | * | | | | | * |
| Creatinine | | * | * | * | | | | | * |
| Cystathionine | * | | * | * | | | | | * |
| Cystine | * (1.25) | | * | * | | * | * | | * |
| Cysteic acid | | | | | | | | * | |
| Ethanolamine | | * | * | * | | | | | * |
| Glutamic acid | * | | * | * | | * | * | * | * |
| Glycine | * | | * | * | | * | * | * | * |
| Histidine | | * | * | * | | * | * | * | * |
| D,L-Homocystine | | * | * | * | | | | | * |
| L,L & allo-Hydroxylysine | | * | * | * | | * | | | * |
| 4-trans-L-Hydroxyproline | * | | * | * | | * (1.25) | | | * |
| Isoleucine | * | | * | * | * | * | * | * | * |
| Leucine | * | | * | * | * | * | * | * | * |
| Lysine | | * | * | * | | * | * | * | * |
| Methionine | * | | * | * | * | * | * | | * |
| Methionine-D,L-sulfoxide | | | | | | * | | | |
| Methionine-D,L-sulfone | | | | | | | | * | |
| 1-Methyl-histidine | | * | * | * | | | | | * |
| 3-Methyl-histidine | | * | * | * | | | | | * |
| Norleucine | | | * | | | | | | |
| Ornithine | | * | * | * | | | | | * |
| Phenylalanine | * | | * | * | * | * | * | * | * |
| Phosphoethanolamine | * | | * | * | | | | | * |
| Phosphoserine | * | | * | * | | | | | * |
| Proline | * | | * | * | | * (1.25) | * | * | * |
| Sarcosine | * | | * | * | | | | | * |
| Serine | * | | * | * | | * | * | * | * |
| Taurine | * | | * | * | | | | | * |
| Threonine | * | | * | * | | * | * | * | * |
| Tryptophan | | * | * | * | | | * | | * |
| Tyrosine | * | | * | * | * | * | * | | * |
| Urea | * | | | * | | | | | * |
| Valine | * | | * | * | | * | * | * | * |

NOTE: Concentration of all the constituents in the Amino Acid standards is 0.25 µmole/mL unless otherwise specified.
 *Concentration of all the constituents is 2.5 µmole/mL unless otherwise specified