# www. LCTech.de

Matrix of the Month



August 2015:

Aflatoxins B/G and **Ochratoxin A in Raisins: Fully Automated with FREESTYLE ThermELUTE™** 




Do you have a special matrix that we should test for mycotoxins Please let us know and write an e-mail to info@LCTech.de!

# Fully Automated Sample Preparation and Analysis

Automate the sample preparation in your laboratory and let the robotic system FREESTYLE ThermELUTE™ handle your samples from raw extract to chromatogram: 24 hours / day for 7 days / week.

Due to the unique technology of thermal denaturation combined with the miniaturisation of the process you need less sample. At the same time you have less solvent consumption and achieve more sensitive and reproducible results in the lower ppt-range and a higher throughput.

# Automated Processing with FREESTYLE ThermELUTE™



food, the analytes can be reliably mesured. Due to the miniaturisation of the process using the SMART immunoaffinity columns

Using the FREESTLYE ThermELUTE™ the results become more sensitive thus even in ranges with very low limits, e.g. for baby

AflaCLEAN SMART and OtaCLEAN SMART less sample, less solvent and less processing time is required.

Raisins, 10 ppb aflatoxins

Black line:

FREESTYLE ThermELUTE™, 2.5 mL sample

Red line:

FREESTYLE ThermELUTE™, 1 mL sample

Orange line:

Manual processing with AflaCLEAN column; 50 mL sample



Extract 20 g of homogenised raisins and 2 g NaCl with 100 mL of the extraction solution (methanol/water, 80/20, v/v) and 50 mL n-hexane for at least 10 minutes for defatting. Filtrate the raw extract.



You can centrifuge it for 10 minutes at 2000 x g facilitating the phase separation.

Aflatoxin: Dilute 7 mL of the filtrate with 43 mL PBS. In case of precipitation or turbidity you can remove them by filtration. Apply max. 10 mL thereof onto the immunoaffinity column AflaCLEAN SMART by FREESTYLE ThermELUTE™.

Ochratoxin: Dilute 10 mL of the filtrate with 40 mL PBS. In case of precipitation or turbidity you can remove them by filtration. Apply max. 10 mL thereof onto the immunoaffinity column OtaCLEAN SMART by FREESTYLE ThermELUTE™. (Shown in the subsequent chromatogram: 3 mL of the filtrate are diluted with 12 mL PBS and 10 mL thereof are applied onto the immunoaffinity column OtaCLEAN SMART by FREESTYLE ThermELUTE™. This represents 0.08 g matrix.

The respective column is automatically loaded, washed with 2 mL water (Aflatoxin: 3 mL/min., Ochratoxin: 1.5 mL/min) and heated by the FREESTYLE system. Afterwards the eluate is quantitatively injected into the sample loop and analysed via HPLC-FLD.

Aflatoxins R/G

### **HPLC Conditions**

	Aliatoxilis b/o	Ochratoxiii A		
HPLC:	Isocratic	Isocratic		
Column oven:	36 °C	40 °C		
Separation column:	RP C18	RP EC 125/3 nucleosil 120-3 C18		
Flow rate:	1.2 mL/min	0.6 mL/min		
Eluent:	HPLC-water/methanol/ acetonitrile (60/30/15 (v/v/v))	HPLC-water/methanol/ acetonitrile (40/55/5 (v/v/v)) + 1 % acetic acid		
Fluorescence detection:	Photochemical derivatisation with UVE	Without derivatisation		
Excitation wavelength:	365 nm	335 nm		
Emission wavelength:	460 nm	465 nm		

Ochratovin A

# Recovery Rates


		Aflatoxin			Ochratoxin A
	B1	B2	G1	G2	ОТА
Standard*	100	100	100	100	100
Recovery rate** Raisins 10 ppb	94	96	94	91	90

SMART columns with tip

**Chromatograms** 




## Chromatograms






Fully automated, sensitive,


Ota standard, 10 ppb (0.8 ng / injection)



Raisins, 10 ppb (0.08 g loaded, spiked with 10 ppb OTA)



Raisins, 1 ppb OTA (0.08 g matrix were injected)



Calibration curve OTA (n=3): (0-1.2 ng/10 mL represents 0-12.5 ppb); LOQ 0.03 ppb (signal/noise 5:1)

These LCTech products were used:

AflaCLEAN SMART immunnoaffinity column for aflatoxins B1, B2, G1 and G2

P/N 12862 / 12863

OtaCLEAN SMART immunnoaffinity column for ochratoxin A

P/N 13346 / 13351

FREESTYLE ThermELUTE™
Robotic system
for sample preparation and
analysis

P/N 12663 / 12668 / 13691

Do you have further questions? Please just write an e-mail to info@LCTech.de!

www. LCTech.de