



Aflatoxins B/G und Ochratoxin A in Tofu

manual and automated



Immunoaffinity columns for aflatoxins B1, B2, G1 and G2

AflaCLEAN and OtaCLEAN are suitable for the analysis and clean-up of various matrices, whether feed or food, from A like amaranth to Z like cinnamon, all matrices can be analysed with reproducible results and highly sensitive to the mycotoxins aflatoxin B/G and/or ochratoxin A. Both the AflaCLEAN and the OtaCLEAN columns are suitable for manual and automated processing, e.g. with the LCTech FREESTYLE SPE robotic system. Simply place the two columns on top of each other and they are ready for the purification of

Processing protocol

Extract 20 g tofu with 2 g sodium chloride in 100 mL methanol/water (80/20 (v/v)) for 30 minutes. Add 50 mL n-hexane during the 30-minute extraction to extract fats and oils. Filter the crude extract, if there is no clear phase separation after filtration, centrifugation at 3000 xg for 10 minutes is helpful. Dilute 2 mL of the methanolic (lower) phase with 12 mL of PBS buffer containing 8 % Tween20, this reduces non-specific binding of matrix components to the IAC column. Load 14 mL onto an immunoaffinity column (AflaCLEAN, or OtaCLEAN), this corresponds to a matrix amount of

Tofu

Tofu is made from white soybeans processed into soy milk. The soy curd produced by denaturing and coagulating proteins is dehydrated and then pressed into blocks. In cereals and legumes such as soybeans, mycotoxins can be produced by incorrect storage conditions and during the drying process. As these can be harmful to humans and animals, soybeans and their products are regularly tested for toxins. In the case of imports, EU border controls show that European limit values are not always complied with. For example, aflatoxin and ochratoxin A levels that are too high repeatedly lead to food being rejected.

Best sample purification for food and feed

Immunoaffinity columns with SMART-columns made by LCTech



aflatoxins B1, B2, G1, G2 and ochratoxin A.

Advantages of the AflaCLEAN column at a glance:

- 1 mL and 3 mL format for homogeneous running speed and optimal clean-up
- Outstanding shelf life up to 24 months at room temperature
- Loading capacity: 150 ng aflatoxin B1 (with OtaCLEAN even: 200 ng ochratoxin A)
- Recoveries: B1 > 90 %, B2 > 80 %, G1 > 90 %, G2 > 60 %
- Suitable for automated processing

0.4 g. Rinse the sample introduction vessel with 5 mL deionised water and load the rinse solution onto the immunoaffinity column. The column is washed again with 5 mL deionised water. After the wash solution has passed through the column, the column is dried by a stream of air and eluted using 2 mL methanol. The methanol should be allowed to soak into the column bed for at least 5 minutes to allow complete denaturation of the antibodies and thus elution of the toxins. The eluate is diluted to runner ratios or injected in small volumes and analysed. The samples can thus be analysed by HPLC fluorescence or by LC-MS/MS.



Conclusion

Tofu, as a soy-based food, can also be contaminated by mycotoxins found in soy. Checking food for residues and contaminants using the immunoaffinity columns AflaCLEAN and OtaCLEAN allows the simple detection and purification of mycotoxins (aflatoxins B/G and ochratoxin A), which as storage fungi are significantly responsible for health damage.

Thus, food and feed can be assessed as safe and of high quality. Even processed and complex foodstuffs can be tested by means of the immunoaffinity columns AflaCLEAN and OtaCLEAN for mycotoxins.

The samples are compatible with HPLC fluorescence or LC-MS/MS analysis. High matrix depletion and concentration of the analytes ensure best results and sensitive analytics according to the strict limits set.



These LCTech products were used:

10514	AflaCLEAN
10515	OtaCLEAN
10519	UVE
10522	HPLC column for mycotoxins
10750	Pre-column holder
10523	Pre-column (Guard)

Recovery rates** Aflatoxins

Aflatoxin	B1	B2	G1	G2
Standard*	100	100	100	100
Tofu AflaCLEAN 10 ppb	93	82	103	98

Recovery rates** Ochratoxin

Ochratoxin A	
Standard*	100
Tofu OtaCLEAN 10 ppb	91

* Standard was set = 100% ,

** Corrected with non-spiked sample / The results are in accordance with the performance specifications of EC 401 / 2006 (section 4.3.1).

Conditions

	Aflatoxins	Ochratoxin A
HPLC	Isocratic	Isocratic
Column oven	36 °C	40 °C
Separation column	RP C18 PN 10522	RP EC 125/3 nucleosil 120-3 C18
Flowrate, SolventI	1.2 mL/min	0.6 mL/min; HPLC-water/methanol/acetonitrile (40/55/5 (v/v/v)+1% acetic acid)
Derivatisation	Photochemical by UVE	without derivatisation
Fluorescence detection		
Excitation wavelength	365 nm	335 nm
Emission wavelength	460 nm	465 nm

Do you have a special request as to which matrix we should test for you?
Contact us by e-mail at: info@LCTech.de