

# Analysis of *PFAS from Soil* Using EluCLEAN PFAS SPE Columns

Authors

Dr. Suman Kharel, LCTech GmbH

Dr. Thomas Gersthagen, LCTech GmbH





# Content

Key Features.....	3
LC Tech Products .....	3
Other Relevant LC Tech Application Notes and Product Information .....	3
1. Introduction.....	4
2. Experimental .....	5
2.1 Sample Preparation.....	5
2.1.1 Sample Extraction.....	5
2.1.2 Solid Phase Extraction.....	5
2.1.3 Evaporation/Concentration .....	6
2.2 Instrumentation.....	6
2.2.1 MS Conditions.....	6
2.2.2 LC Instrument Conditions .....	6
3. Results.....	7
3.1 Analyte, their Retention Times and Corresponding Isotope Dilution Standard.....	9
3.2 Recovery Rates and RSD % of 40 PFAS.....	11
4. Conclusion.....	13
5. References.....	13



## Key Features

- Excellent recovery rates and low standard deviations for 40 PFAS analytes according to US EPA 1633 (2nd draft)
- No detectable PFAS background contamination
- Only one SPE column for clean-up and enrichment needed
- EluCLEAN PFAS – WAX/GCB SPE column can be used equivalently to the SPE cartridge + dispersive graphitized carbon black used in US EPA 1633 (2nd draft)
- EluCLEAN PFAS – WAX/GCB SPE column contains a weak anion exchanger, mixed-mode polymeric sorbent with an pKa above 8 with optimised parameters for PFAS enrichment, suitable for use in US EPA 1633 (2nd draft) and DIN 38414-14
- EluCLEAN PFAS – WAX/GCB SPE column with higher sorbent amount for samples with higher matrix load is also suitable for DoD/QSM 5.1/5.3 and DIN 38414-14

## LCTech Products

### SPE cartridges

Part No.: 20821, 20822, 20823

EluCLEAN PFAS – WAX/GCB

Sorbent 1: Weak Anion Exchanger, Mixed-Mode Polymeric Sorbent (WAX),

Sorbent 2: Graphitized Carbon Black (GCB)

150/10 mg/6 mL

Part No.: 20831, 20832, 20833

EluCLEAN PFAS – WAX/GCB

Sorbent1: Weak Anion Exchanger, Mixed-Mode Polymeric Sorbent (WAX)

Sorbent 2: Graphitized Carbon Black (GCB)

200/50 mg/6 mL

## Other Relevant LCTech Application Notes and Product Information

AN0052 Analysis of PFAS from Drinking Water Using EluCLEAN PFAS - SPE columns

AN0054 Analysis of PFAS from Drinking Water Using Automated FREESTYLE XANA-PFAS System and EluCLEAN PFAS - SPE Columns

AN0045 D-EVA – Automated EVApOration of PFAS compliant to US-EPA 537.1

Flyer MIX-TRACTION (P/N 20444)



# 1. Introduction

Per- and polyfluorinated alkyl substances (PFAS) products have been in use for more than 60 years. They get into the environment during their manufacturing process, usages and disposal. Research has revealed the high toxicity of PFAS compounds and thus the resulting need to regulate the substances. Therefore the analytical interest in these compounds has rapidly increased in the last few years. The current and upcoming regulations in the EU and US make it necessary to test soil from various locations for PFAS content. Different methods for PFAS analysis in the EU and US exist. For example, DIN 38414-14, US-EPA 1633 (2nd draft) [1] and DoD/QSM 5.1/5.3. All methods require solid phase extraction (SPE) prior to liquid chromatography-tandem mass spectrometry (LC-MS/MS) analysis. All methods apply SPE cartridges containing a weak anion exchanger, mixed-mode polymeric sorbent, whereas the US EPA 1633 (2nd draft) and the DoD/QSM 5.1/5.3 additionally use a dispersive clean-up step depending on the type of matrices.

In this application note a new single SPE cartridge solution with a PFAS enrichment optimised polymeric sorbent is presented. The EluCLEAN PFAS – WAX/GCB SPE column contains 150 mg of a weak anion exchanger, mixed-mode polymeric sorbent mixed with 10 mg of graphitized carbon black. The SPE cartridge shows excellent recovery rates in combination with low standard deviations and is therefore ideally suited for SPE of PFAS from soil and other environmental matrices. It can equivalently replace the dispersive clean-up step + WAX SPE cartridge used in US EPA 1633 (2nd draft). The single cartridge solution saves time and costs in PFAS analysis.




Figure 1. EluCLEAN PFAS – SPE columns



## 2. Experimental

### 2.1 Sample Preparation

#### 2.1.1 Sample Extraction

An aliquot of soil (5 g dry weight) was filled into a 50 mL polypropylene centrifuge tube. 10 mL of 0.3 % methanolic ammonium hydroxide were added to each centrifuge tube. The sample was vortexed and then shaken for 5 minutes at 800 rpm with the MIX-TRACTION system. After centrifugation at 2800 rpm for 10 minutes the supernatant was transferred into a clean 50 mL polypropylene centrifuge tube. The process was repeated two more times with 15 mL and 5 mL 0.3 % methanolic ammonium hydroxide solution.

The samples were then diluted to 35 mL by adding 5 mL reagent water. After concentration and removal of methanol with the D-EVA Rotational Vacuum Concentrator, the concentrated final volume was approx. 7 mL per sample. The samples were then filled with reagent water to get a final volume of 45 mL. After vortexing the pH was measured and adjusted to  $6.5 \pm 0.5$  with 50 % formic acid and/or 30 % aqueous ammonium hydroxide as required.

The samples were spiked with 40 PFAS (PFAC-MXF-J, Wellington Lab) and 24 mass labelled extraction standard (PFAC-HIF-ES, Wellington Lab). The resulting final concentration were 0.5 - 40  $\mu\text{g}/\text{kg}$  soil for both standards.

#### 2.1.2 Solid Phase Extraction

EluCLEAN PFAS – WAX/GCB (150 / 10 mg / 6 mL) SPE cartridges were placed on an EluVac vacuum manifold and 25 mL reservoirs were placed on top of each cartridge. The cartridges were pre-conditioned with 15 mL of 1 % methanolic ammonium hydroxide followed by 5 mL of 0.3 M formic acid (cartridges were not allowed to run dry). The samples were applied with approx. 5 mL/min (using vacuum). The cartridges and the reservoirs were washed with 5 mL reagent water (twice) followed by 5 mL of 1:1 0.1 M formic acid/methanol using vacuum. The cartridges were then dried with high vacuum for 15 sec. The 50 mL sample polypropylene centrifuge tube was rinsed using 5 mL of 1 % methanolic ammonium hydroxide then, using a glass pipette. The rinsate was transferred to the reservoir, washing the walls of the reservoir. Vacuum was slowly applied to pull the elution solvent through the cartridge into the 15 mL polypropylene centrifuge tubes.

25  $\mu\text{L}$  of concentrated acetic acid and 10  $\mu\text{L}$  NIS (MPFAC-HIF-IS) solution were added to each sample eluate and vortexed. The samples were filtered through a syringe filter (25 mm filter, 0.2  $\mu\text{m}$  nylon membrane) into a clean 15 mL polypropylene centrifuge tube.



### 2.1.3 Evaporation/Concentration

All samples were evaporated to 500  $\mu$ L - 1 mL using D-EVA Rotational Vacuum Concentrator (temperature: 45°C, vacuum: 20 mbar) and transferred into a 1.5 mL polypropylene vial and kept at 0 – 4 °C for LC-MS/MS analysis.

## 2.2 Instrumentation

### 2.2.1 MS Conditions

Table 1. MS Conditions

Parameter	Value
MS	TSQ Quantis (Thermo)
Polarity	Negative
Spray voltage	2500 V
Sheath Gas	50 Arb
Aux Gas	10
CID Gas	2 mTorr
Ion transfer tube temp	325 °C
Vaporizer Temp	300 °C
Q1 resolution	0.7 FWHM
Q3 resolution	1.2 FWHM
Cycle time	0.5 sec
Chromatographic peak width	6 sec

### 2.2.2 LC Instrument Conditions

Table 2. LC Conditions

Parameter	Value	
LC	Thermo Scientific Vanquish Flex UHPLC system	
Analytical column	Accucore RP-MS, 2.1*100 mM, 2.6 $\mu$ m	
Delay column	Agilent ZOBAX Eclipse plus C18, 4.6*50 mm; 3.5 $\mu$ m	
Column temperature	45 °C	
Injection volume	5 $\mu$ L	
Mobile Phase	A) 20 mM ammonium acetate H <sub>2</sub> O with 2 % MeOH and 0.1 % acetic acid B) 20 mM ammonium acetate MeOH with 2 % H <sub>2</sub> O and 0.1 % acetic acid	
Gradient Flow rate	0.5 mL/min	
	Time (min)	%
	0	0
	1	30
	6	45
Gradient	13	80
	14	95
	17	95
	20	95
	22	0
	25	0



### 3. Results




Figure 2. Chromatogram of 40 PFAS, 24 labelled extraction standards and 7 internal standards after extraction with EluCLEAN PFAS- WAX/GCB SPE column from soil spiked samples (spiked concentration = 0.5-40 µg/kg soil)




Figure 2. Chromatogram of 40 PFAS, 24 labelled extraction standards and 7 internal standards after extraction with EluCLEAR PFAS-WAX/GCB SPE column from soil spiked samples (spiked concentration = 0.5-40 µg/kg soil)





### 3.1 Analyte, their Retention Times and Corresponding Isotope Dilution Standard

Table 3. Overview retention times and corresponding isotope dilution standard

Analyte	Retention time (min)	Isotope Dilution Standard
PFBA	2.80	MPFBA
PFMPA	3.33	MPFBA
PFPeA	4.51	MPFPeA
PFMBA	5.21	MPFBS
L-PFBS	5.19	MPFBS
PFEESA	6.14	MPFBS
NFDHA	6.55	MPFBS
4:2FTS	6.66	M2-4:2FTS
PFHxA	6.98	M5PFHxA
L-PFPeS	7.59	M5PFHxA
3:3 FTCA	7.93	M5PFHxA
HPFO-DA	7.78	M3HPFO-DA
PFHpA	9.30	M4PFHpA
PFHxSK	9.61	M3PFHxSK
ADONA	9.56	M3PFHxSK
6:2FTS	10.71	M2-6:2FTS
PFOA	10.82	M8PFOA
L-PFHpS	10.97	M9PFNA
5:3 FTCA	11.74	M9PFNA
PFNA	11.96	M9PFNA
PFOSK	12.03	M8PFOS
9Cl-PF3ONS	12.51	M9PFNA
8:2FTS	12.81	M2-8:2FTS
PFDA	12.86	M6PFDA
L-PFNS	12.89	M7PFUnA
PFUnA	13.63	M7PFUnA
7:3 FTCA	13.73	M8FOSA
L-PFDS	13.63	M8FOSA
N-MeFOSAA	13.60	d3-N-MeFOSAA
FOSA	13.77	M8FOSA
11Cl-PF3OUdS	13.98	M8FOSA
N-EtFOSAA	14.04	d5-N-EtFOSAA
PFDoA	14.28	MPFDoA
N-MeFOSA	14.67	d-N-MeFOSA
N-MeFOSE	14.60	d7-N-MeFOSE



L-PFDoS	14.77	MPFDoA
PFTrDA	14.82	M2PFDA
N-EtFOSA	15.05	d-N-EtFOSA
N-EtFOSE	14.97	d9-N-EtFOSE
PFTeDA	15.15	MPFTeDA



### 3.2 Recovery Rates and RSD % of 40 PFAS

Table 4. Recovery rates and RSD

Analytes	Recovery rate (%)	RSD (%)
PFBA	99	9
PFMPA	102	11
PFPeA	96	9
PFMBA	101	11
L-PFBS	93	15
PFEESA	98	11
NFDHA	96	12
4:2FTS	104	13
PFHxA	103	12
L-PFPeS	94	9
3:3 FTCA	91	5
HPFO-DA	86	8
PFHpA	107	9
PFHxSK	94	14
ADONA	93	14
6:2FTS	108	11
PFOA	108	3
L-PFHpS	89	18
5:3 FTCA	99	11
PFNA	106	9
PFOSK	88	10
9Cl-PF3ONS	96	6
8:2FTS	93	8
PFDA	101	11
L-PFNS	111	12
PFUnA	103	13
7:3 FTCA	115	7
L-PFDS	104	20
N-MeFOSAA	92	11
FOSA	94	11
11Cl-PF3OUdS	89	5



N-EtFOSAA	98	8
PFDaA	111	9
N-MeFOSA	92	13
N-MeFOSE	96	17
L-PFDoS	110	13
PFTrDA	113	14
N-EtFOSA	106	18
N-EtFOSE	95	17
PFTeDA	93	13




Figure 3. Recovery rates of 40 PFAS (listed in US EPA method 1633 2<sup>nd</sup> draft) from soil matrix extracted with EluCLEAR PFAS- WAX/GCB SPE column (n = 4, spiked concentration = 0.5-40 µg/kg soil )



## 4. Conclusion

EluCLEAN PFAS – WAX/GCB cartridges have no detectable PFAS background contamination. The cartridge EluCLEAN PFAS – WAX/GCB 150/10 mg is fulfilling the required performance of US EPA 1633 (2nd draft). Recoveries for 2 - 40 µg/kg soil samples with a range of 88 - 115 % were very well in between the acceptable criteria of 70 - 130 %. Therefore, the desired accuracy is given. Precision, measured by % RSD of replicate extracts, was also well within the range of requirements, with all RSDs below 20 %.

EluCLEAN PFAS – WAX/GCB cartridges are therefore ideally suited to be used for the enrichment and clean-up of PFAS from soil and other environmental matrices. It can replace the dispersive carbon clean-up + WAX cartridge or other two cartridge solutions with only one column. This is the time and cost saving alternative for every laboratory.

## 5. References

[1] 2nd Draft Method 1633: Analysis of Per- and Polyfluoroalkyl Substances (PFAS) in Aqueous, Solid, Biosolids, and Tissue Samples by LC-MS/MS; EPA Document No. EPA 821-D-22-001, June 2022.

Any Questions?  
Do not hesitate to contact us: