™

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Dear Sir or Madam,

With your decision to buy the vacuum manifold EluVac you have made an excellent choice. Congratulations!

The concept of this system is the consistent result of experiences that diverse laboratories have made with solid phase extraction over the years.

Using the EluVac both small and large volume sample are easily and quickly handled with the same system.

Moreover eluate vials and racks needn't be put into the vacuum system during the process. Simply turn the lid and all positions are changed from "WASTE" to "COLLECT" position.

Much attention was paid to high quality materials and components, e.g. for a precise vacuum regulation. The whole system is designed for daily and at the same time long-term use.

With appropriate maintenance and careful handling you will be able to use the vacuum manifold over a long period of time without service.

In case you should have questions or special requirements, e.g. a special rack, the LCTech team and your distributor, respectively, will be pleased to support you.

The LCTech-Team wishes you all the best with your new EluVac!

P.S. Do you have changes or remarks to this manual? Please directly contact us with an e-mail to service@lctech.de.

1 Intended Purpose

The SPE vacuum manifold EluVac is designed for handling SPE and immunoaffinity columns in chemical laboratories with a vacuum to be regulated by the user.

The used SPE columns must have a male Luer connection at the column outlet.
Up to 20 columns can be handled in parallel.

Depending on the expected waste volume the EluVac can be used with two different versions for drainage.

2 Safety Instructions

Please mind these safety instructions!

Please read the manual and the safety instructions.

The system may only be operated by trained, competent staff in appropriate laboratories. Please pay attention to common basic rules for careful work in laboratories.

During operation appropriate protective clothing (e.g. safety glasses, gloves, laboratory coat and shoes) have to be put on. Any physical contact with chemicals has to be avoided.

Take care of an appropriate ventilation.

The system or parts of it must not be used with persons or animals.

Check that the system is complete and free of defects before operation.

Damaged system parts (cracks, corrosion, etc.) have to be exchanged before the system is put into operation.

Please familiarize with vacuum systems.

WARNING: Implosion hazard. The vacuum should not exceed 22 in. Hg (~ 250 mbar absolute).

The vacuum must not exceed the defined limit, as an implosion hazard might occur.

If vacuum pumps are used, their specifications have to be suitable concerning chemical compatibility and vacuum regulation.

Please comply with the advices for connecting and assembling the system.

Be sure that the individual components are assembled according to the manual and that exclusively appropriate materials and vials are used.

The EluVac system was designed for room temperature between +10 °C and +35 °C (50 °F to 95 °F) during operation.

Compounds that are highly oxidizing, explosive, corrosive or otherwise instable, must not be used in the vacuum system.

All liquid and solid wastes have to be collected and disposed of in an appropriate way.

If a Woulff bottle or a similar external waste container or the 600 mL beaker in the system is used, be sure that it is emptied before operation and that the applied volume does not exceed the collection volume. Moreover be sure that the Woulff bottle or other used glassware is appropriately coated (implosion prevention!).


If low-boiling solvents are used, solvent vapour may be produced. Be sure it does not discharge or inflame. Be sure that no open source of ignition is near the system.

Clean the system from time to time and store it properly and safely when it is not used.

Never dispose of laboratory systems or dangerous materials. The materials have to be recycled according to valid regulations or disposed of in an environment-friendly way.

WARNING! Never ignore these safety instructions, as this may result in serious injuries and damages.


3 Basic System


For very frequent applications special sets are available containing all necessary components (see Ordering Information).

4 Assembling the Racks


Carefully take out the single components of the packaging and assemble them as shown here:

<p>Rack for eluate collection in 4 mL vials, up to 20 positions</p>	<p>Part N° 11105</p>
 <p>Upper bar (11102) and set screw (10464)</p> <p>Upper disc (11099)</p> <p>Middle disc (11099)*</p> <p>Middle bar (11103)</p> <p>Lower disc (11100)*</p> <p>Lower bar (11104) and set screw (10464)</p> <p>* between middle and lower disc two separation rings are added (11246).</p>	
<p>Rack for eluate collection in 60 mL vials, up to 10 positions</p>	<p>Part N° 11118</p>


Rack for eluate collection in rotary evaporator bottles,
up to 6 positions

Part N° 11112


Rack for eluate collection in 100 mL vials,
up to 10 positions

Part N° 12767

Rack for eluate collection in 25 mL round bottom flasks,
up to 10 positions

Part N° 13465


5 Assembling the Holders for Sample Reservoirs

The holders are the ideal completion if you want to use large volume sample reservoirs.


Holder for sample reservoirs, 10 positions	Part N° 11098
 <ul style="list-style-type: none"> <li data-bbox="938 515 1364 548">Screw (10722) with washer (10411) <li data-bbox="938 560 1165 593">Upper disc (11094) <li data-bbox="938 772 1077 806">Bar (11095) <li data-bbox="938 985 1165 1019">Lower disc (11096) <li data-bbox="938 1064 1364 1097">Screw (10813) with washer (10411) 	
Holder for sample reservoirs, 6 positions	Part N° 11560
 <ul style="list-style-type: none"> <li data-bbox="938 1276 1364 1310">Screw (10722) with washer (10411) <li data-bbox="938 1321 1173 1355">Upper disc (11570) <li data-bbox="938 1444 1085 1478">Bar (11095) <li data-bbox="938 1635 1396 1668">Screw nut (10470) and screw (10722) <li data-bbox="938 1736 1173 1769">Lower disc (11571) <li data-bbox="223 1825 319 1892">Screw (10813) <li data-bbox="938 1881 1117 1915">Screw (10808) 	

6 Operating the Manifold

- Be sure that waste fluids are collected. Use one of two possibilities for it:
 - Pull the central drain tube into the central boring and connect it to a tubing and an external waste container (e.g. a 10 L Woulff bottle)


- Alternatively: Place a waste container (e.g. the beaker delivered with the system) centrally into the EluVac

- Put the rack for eluate collection into the basic system according to the borings
- Put the eluate vials into the rack for eluate collection
- Close the basic system with the inner lid
- Place the transfer tips into the positions to be used
- Put up the upper lid according to the borings

- Please note the labels „WASTE“ and „COLLECT“, for a switch slightly lift the lid and turn it into the right direction


- Place the column valves at those positions in the upper lid to which eluate vials are placed below
- Close the unused positions with the port sealing plugs
- If you use Luer adapters: Put them onto your columns

- If you use a holder for sample reservoirs: Put it upon the system acc. to the borings
- Apply the vacuum: Switch on the vacuum pump and push the lid down equally at both handle bars
- Put up the samples
- Regulation of the vacuum: Coarse regulation is done with the vacuum valve, fine regulation is done with the column valves. Each column can be regulated individually

7 Set-Up for Large Volume Samples

For large sample volumes this set-up is recommended:

- Several samples can be handled in parallel. The individual samples can have up to several litres.
- Via tubing the sample is aspirated from the sample bottle and applied to the column, which is connected with an adapter.
- The waste fluid is drained off through a central drain tube and the manifold ...
- ... and collected in an external, large container, e.g. a 10 L Woulff bottle.
- The vacuum pump generates the vacuum.

This set-up is according to the EluVac set with the part N° 11739 (details see Ordering Information).

8 Ordering Information

BASIC SYSTEM		Part N°
EluVac Vacuum Manifold (W x D x H: 30 x 30 x 21 cm) incl. vacuum regulation unit, beaker / drain tube with tubing, 20 column valves, 20 port sealing plugs		11048
COMPLETE-SETS (for common applications, without vacuum pump or consumables)		
for large sample volumes		11739
for the determination of the hydrocarbon oil index		11146
for the analysis of mycotoxins		12415
for pear shaped rotary evaporator bottles, NS 29/32 (up to 100 mL)		11147
RACKS FOR ELUATE COLLECTION		
for 4 mL vials (V0004), 20 positions		11105
for 60 mL vials (F060), 10 positions		11118
for 120 mL vials (F100), 10 positions		12767
for 100 mL rotary evaporator bottles, NS 29/32 (11106), 6 positions		11112
HOLDER FOR SAMPLE RESERVOIRS		
10 positions		11098
6 positions, only for rack with rotary evaporator bottles (11112)		11560

ACCESSORIES

Sample reservoir (2 pieces / pack) for 100 mL sample, made of glass, reusable and dishwasher proof, suitable for LCTech Florisil columns, incl. screw cap, PTFE seals	9370-1
Sample reservoir (2 pieces / pack) for 100 mL sample, made of glass, reusable and dishwasher proof, suitable for LCTech immunoaffinity columns, incl. screw cap and PTFE seals	10896
Seals for sample reservoir (9370-1), PTFE (10 pieces / pack)	9373-2
Seals for sample reservoir (10896), PTFE (10 pieces / pack)	10909
Woulff bottle, 10 L, incl. special screw cap and 2 m tubing	F10000B
Adapter for SPE columns, for sucking larger sample volumes of an external container, suitable for all 1, 3, 6 or 15 mL columns, incl. 1.5 m tubing	11734
Adapter from ELUFIX column to Luer connection (10 pcs./pack)	11043
Transfer tips (yellow, made for EluVac) (50 pcs./pack)	11597

VACUUM PUMPS

ME 2C, one stage, "Chemistry", incl. cable, plug and instructions 230 V, 50-60 Hz	11125
100 – 120 / 200 – 230 V, 50-60 Hz	11126
ME 2, one stage, „Aluminum“, incl. cable, plug and instructions 230 V, 50-60 Hz	11123
100 – 120 / 200 – 230 V, 50-60 Hz	11124

SAMPLE VIALS OF GLASS

4 mL vials, 100 pieces / pack	V0004
60 mL vials, 100 pieces / pack	F060
120 mL vials, 2 pieces / pack	F100-D
100 mL rotary evaporator bottles, NS 29/32, pear shaped, 1 piece	11106

Any Questions?
Do not hesitate to contact us:

Limitation of Liability

In no event will LCTech be liable to any part for any indirect, incidental, consequential, exemplary, special or punitive damages, including any loss of profit, revenue, business opportunity or data, arising from or relating to this equipment, whether in contract, in tort, negligence or otherwise, even if LCTech know, should have known or has been advised of the possibility of such damages. LCTech's total cumulative liability arising from or related to this equipment, whether in contract, in tort or otherwise, will not exceed the price actually paid by you for this equipment.