

Iso-Alpha Acids in Beer How About Automating Sample Preparation?

Dear Sir or Madam,

March 15, 2010

Are you concerned with the analysis of iso-alpha acids in beer samples?

We have successfully tested our laboratory system AcceCLEAN for the analysis of iso-alpha acids in beer in cooperation with a well-known company. Maybe the AcceCLEAN could also be a good opportunity for you to automate major parts of this analysis - at a price which might be lower than you expect it to be.

Using the AcceCLEAN **you can separate the iso-alpha acids from the beer sample with solid phase extraction:**

Beer sample



The AcceCLEAN fully automatically handles the samples via C18 SPE columns

Three samples in parallel,
30 in one run



Result:
Iso-alpha acids in the eluate (water/methanol), ready for analysis

Do you wish to gather more information? Please simply contact me!

All the best!



Christine Grossmann, Sales Manager

*P.S. We have developed this application in close cooperation with one of our costumers.
Are you looking for the right way how to automate your SPE application? Please simply ask us!*