



AflaCLEAN M1 SMART, Immunoaffinity Columns

Ordering Numbers: 14246, 14248, 14249

Packaging Sizes: 5, 100, and 1000 columns

The immunoaffinity column (polypropylene) for sample preparation containing a gel material loaded with monoclonal antibodies against the mycotoxin aflatoxin M1 in a storage buffer.

The column may only be used for food testing and in accordance with the LCTech testing procedure (extraction and clean-up) and HPLC with fluorescence detection.

The columns are designed for a one-time use only; defective test results as a result of any other use are not covered by the responsibility of LCTech.

Store the columns at 4 °C and do not freeze.

The shelf life is 9 month; for expiry date see package label.

The minimum recoveries in percent are:

Aflatoxin M1 > 90

The maximum capacity of the column is 100 ng aflatoxin M1.

Date of issue: 17/04/10



The company LCTech GmbH is certified according to ISO 9001:2015

