



Biotin Clean-Up Column

Ordering Numbers: 14099, 14100, 14101

Packaging Sizes: 5, 25, and 500 columns

The affinity column (polypropylene) for sample preparation containing a gel material loaded with specific molecules to bind biotin in a storage buffer.

The column may only be used for food testing and in accordance with testing procedure (extraction and clean-up) and HPLC with UV or fluorescence detection.

The columns are designed for a one-time use only; defective test results as a result of any other use are not covered by the responsibility of LCTech.

Store the columns at 4 °C and do not freeze.

The shelf life is 9 month; for expiry date see package label.

The minimum recoveries in percent are:

Biotin: ≥ 90

The maximum capacity of the column is 1000 ng biotin and can be lower for biotin derivatives.

Date of issue: 17/04/11



The company LCTech GmbH is certified according to ISO 9001:2015

