



CrossTOX®, Clean-up Columns for Multi-Mycotoxin Analysis

For Automated FREESTYLE SPE / QuEChERS Processing

Ordering Numbers: 17901

Packaging Sizes: 100 columns

The clean-up column (polypropylene) for sample preparation contains a solid phase bed material for purification of mycotoxins from various food and feed material.

The column may only be used for food and feed testing and in accordance with the LCTech testing procedure (extraction and clean-up).

The columns are designed for a one-time use only; erroneous test results because of any other use are not covered by the responsibility of LCTech.

Store the columns at ambient temperature and do not freeze. Avoid contact to humidity.



The company LCTech GmbH is certified according to ISO 9001:2015

