



## ZeaCLEAN SMART, Affinity Columns

**Ordering Numbers:** 14740, 14741, 14742

**Packaging Sizes:** 5, 100, and 1000 columns

The affinity column (polypropylene) for sample preparation contains a bed material for purification and concentration of zearalenone.

The column may only be used for food testing and in accordance with the LCTech testing procedure (extraction and clean-up) and HPLC with fluorescence detection.

The columns are designed for a one-time use only; defective test results as a result of any other use are not covered by the responsibility of LCTech.

Store the columns at 8-30 °C and do not freeze.

The shelf life is 12 month; for expiry date see package label.

The minimum recoveries in percent are:

Zearalenone:  $\geq 90$

The maximum capacity of the column is 2000 ng zearalenone.

Date of issue: 17/04/11



The company LCTech GmbH is certified according to ISO 9001:2015

