

Trade fairs cancelled. And now? LCTech presents its product news.

BE PART OF IT



Login 

German  

Print 

Search

*LCTech-en*

*News*

*Videos*

*FREESTYLE Videos*

# FREESTYLE Videos

Multiple Applications on Automated Clean-up System FREESTYLE (SPE, GPC, EVaporation, ...)

One system fits it all - High throughput or apart from routine or standard tasks. FREESTYLE offers solutions tailored to your requirements.

# Fast and Easy SPE-Method Transfer on Automated Clean-up System FREESTYLE

See for yourself how easy the implementation of automated sample preparation in the laboratory can be.

## FREESTYLE XANA

With 24/7 operation, the parallelisation of the processing procedures and the simultaneous processing of 3 samples, a high sample throughput is achieved. While the system is reliably processing the samples, the user of the system has time to focus on other tasks.

## FREESTYLE GPC/EVaporation User Video

User Video about automated sample preparation for residue analysis.

>> [Click here for the video \(Berghof Analytik + Umweltengineering GmbH\)](#)

## FREESTYLE Bidirectional SPE (BD-SPE)

This video shows a new procedure of sample for the automated sample preparation:

Bidirectional SPE. Thus, even very difficult matrices with smallest sample volumes, like brain, can be processed automatically without supervision over 24 hours and 7 days per week.

## FREESTYLE SPE

Many methods, many solvents, many SPE columns - one ingenious principle.

Solid Phase Extraction (SPE) can be simply automated using FREESTYLE SPE.

## **FREESTYLE ThermELUTE™**

The robotic system FREESTYLE ThermELUTE™ provides fully automated mycotoxin analysis: reliable, reproducible and highly sensitive.



## **FREESTYLE EVaporation**

Automated concentration of all types of organic solvents to a defined end volume between 5.0 mL and 0.2 mL; precision even in unobserved moments.

## FREESTYLE GPC

Automated Gel Permeation Chromatography, easy in handling, efficient and applicable for all food and feed samples.

Please contact us

 +49 8082 

2717-0info@LCTech.de

---

Brochure

[General FREESTYLE brochure \(pdf | 3 MB\)](#)

LCTech on YouTube

You can also watch our videos  
on YouTube.

Video of analytica 2016

Have a look at LCTech in video  
for "food analysis at analytica  
2016" on YouTube.

© 2020 LCTech GmbH

[Disclaimer](#) |  
[Data Security Statement](#) |  
[Terms of Service](#) |  
[Sitemap](#) |