

Trade fairs cancelled. And now? LCTech presents its product news.

BE PART OF IT

Login 

German  

Print 

Search

LCTech-en

Pickering

Artificial Body Fluids

Artificial Sebum and Cerumen

Artificial Sebum and Cerumen

Artificial Sebum

Sebum is an oily secretion produced by sebaceous glands, which spreads over the hair and skin for waterproofing purposes. The formulation of artificial sebum provides the reliability, reproducibility and convenience needed for testing.

Manufactured according to ASTM designation D4265-14 or D4265-98

Used for evaluating the efficacy of home laundry products and conditions to remove stains from fabric. Can also be used in combination with [Artificial Eccrine Perspiration](#) to mimic non-exercise induced skin surface film liquids.

According to	AATCC Standard Test Method 130-2010
Shelf life	one year, refrigerated
Condition	ready-to-use solution
Quantity	25 g
Part Number	1700-0700

Artificial Eccrine Perspiration - Artificial Sebum Emulsion

Combination of [artificial eccrine perspiration](#) and artificial sebum to mimic non-exercise induced skin surface film liquids.

Manufactured according to ASTM D4265-98

Used for evaluating the efficacy of home laundry products and conditions to remove stains from fabric. This solution can be used to test any product that will come in contact with human sweat.

According to	inspired by ASTM D4265-98
Shelf life	one year, refrigerated
Condition	ready-to-use solution
Quantity	200 mL
Part Number	1700-0702

Artificial Cerumen



Cerumen, also known as ear wax, is a waxy substance secreted in the ear canal that protects and lubricates the ear canal and assists in cleaning by trapping dirt and dead skin cells.

The artificial ear wax can be used for testing of hearing aids, ear buds and other electronic devices meant to be used in the ears.

Shelf life: one year, refrigerated

Quantity: 50 g (P/N 1700-0701)

200 g (P/N 1700-0711)

Please contact us

[+49 8082 2717-0](tel:+49808227170) info@LCTech.de

Brochures

[Test_Solutions_Catalog.pdf \(pdf | 3 MB \)](#)

Special Website

Please find detailed information about the test solutions on the special designed Pickering website.

© 2021 LCTech GmbH

[Disclaimer](#) |

[Data Security Statement](#) |

[Terms of Service](#) |

[Sitemap](#) |