

Login 

German  

Print 

Search

LCTech-en

Pickering

Artificial Body Fluids

Artificial Sebum and Cerumen

Artificial Sebum and Cerumen

Artificial Sebum

Sebum is an oily secretion produced by sebaceous glands, which spreads over the hair and skin for waterproofing purposes. The formulation of artificial sebum provides the reliability, reproducibility and convenience needed for testing.

Manufactured according to ASTM designation D4265-14 or D4265-98

Used for evaluating the efficacy of home laundry products and conditions to remove stains from fabric. Can also be used in combination with [Artificial Eccrine Perspiration](#) to mimic non-exercise induced skin surface film liquids.

According to AATCC Standard Test Method 130-2010

Shelf life one year, refrigerated

Condition ready-to-use solution

Quantity 25 g

Artificial Eccrine Perspiration - Artificial Sebum Emulsion

Combination of [artificial eccrine perspiration](#) and artificial sebum to mimic non-exercise induced skin surface film liquids.

Manufactured according to ASTM D4265-98

Used for evaluating the efficacy of home laundry products and conditions to remove stains from fabric. This solution can be used to test any product that will come in contact with human sweat.

According to	inspired by ASTM D4265-98
Shelf life	one year, refrigerated
Condition	ready-to-use solution
Quantity	200 mL

Artificial Cerumen



Cerumen, also known as ear wax, is a waxy substance secreted in the ear canal that protects and lubricates the ear canal and assists in cleaning by trapping dirt and dead skin cells.

The artificial ear wax can be **used for** testing of hearing aids, ear buds and other electronic devices meant to be used in the ears.

Quantity: 50 g, 200 g

Shelf life: one year, refrigerated

Please contact us

+49 8082  

2717-0info@LCTech.de

Brochures

[Test Solutions Catalog.pdf \(pdf | 1.6 MB\)](#)

Special Website

Please find detailed information about the test solutions on the special designed Pickering website.

© 2019 LCTech GmbH

[Disclaimer](#)
[Data Security Statement](#)
[Terms of Service](#)
[Sitemap](#)

|
|
|
|