

Trade fairs cancelled. And now? LCTech presents its product news.

BE PART OF IT

Login 

German  

Print 

Search

*LCTech-en*

*Pickering*

*Artificial Body Fluids*

*Artificial Urine*

# Artificial Urine

Artificial Urine for Testing Sterile Urethral Catheters



According to	DIN EN 1616:1999
pH value	6.6
Storage condition	frozen
Condition	ready-to-use solution, not stabilized
Quantity	200 mL

Part number 1700-0017

### Artificial Urine for Corrosion Testing of Urological Implants

**Designed for** testing metallic biomaterials used to produce urological implants and catheters.

pH value	6.0
Storage condition	room temperature
Condition	ready-to-use solution, stabilized
Quantity	200 mL
Part number	1700-0016

### Artificial Urine Medium for Growing Urological Pathogens

This formulation supports growth of wide range of urinary pathogens and it is also capable of forming crystals similar to these found in natural urinary tract infections.

**Used for** clinical studies and product testing such as negative controls in laboratory testing.

pH value	6.5
Storage condition	room temperature
Condition	ready-to-use solution
Quantity	200 mL

The artificial urine is available in a stabilized version (1700-0600) and non-stabilized version (1700-0018).

Please contact us

 +49 8082 2717-0  [info@LCTech.de](mailto:info@LCTech.de)

---

Brochures

[Test\\_Solutions\\_Catalog.pdf \( pdf | 3 MB \)](#)

Special Website

Please find detailed information about the test solutions on the special designed Pickering website.

© 2021 LCTech GmbH

[Disclaimer](#) |

[Data Security Statement](#) |

[Terms of Service](#) |

[Sitemap](#) |