

Trade fairs cancelled. And now? LCTech presents its product news.

BE PART OF IT

Login 

German  

Print 

Search

LCTech-en

Products

EluVac Vacuum Manifold

Sample Reservoirs

Sample Reservoirs

The 100 mL sample reservoirs are used for the collection of the sample and are suitable for all LCTech stands like e.g. the [vacuum manifold EluVac](#).

The sample reservoirs are suitable for all kinds of 3 mL [SPE columns](#). They can be combined with both the LCTech [immunoaffinity columns](#) and the [Florisil columns](#).

They are made of DURAN glass, lab washer-proof and reusable. The reservoirs are delivered with screwing's and PTFE seals.



LCTech sample reservoirs for large sample volumes - combined with a immunoaffinity column (left) and with a Florisil column on the vacuum manifold EluVac (right).

Please contact us

Brochures

[Mycotoxin brochure \(pdf | 3 MB \)](#)

EluVac racks

[Overview of racks for eluate collection](#)

© 2021 LCTech GmbH

[Disclaimer](#) |

[Data Security Statement](#) |

[Terms of Service](#) |

[Sitemap](#) |