

Trade fairs cancelled. And now? LCTech presents its product news.

BE PART OF IT

Login 

German  

Print 

Search

LCTech-en

Products

GPC Columns

GPC Columns

Functionality Through and Through



The LCTech GPC columns are designed for easy and self-explanatory handling. Due to their robustness they can be used for every-day analysis in laboratories. The columns are available filled (with Bio-Beads) or unfilled as well as in different lengths and filling heights. They can be used with a variety of solvents.

To re-fill the columns directly on site, LCTech offers suitable kits for a simple and fast procedure. The columns can be easily taken apart, cleaned and put together again. Owing to their refined construction, the columns will be instantly tight and ready for use.



Easy refilling of the columns

Further advantages at a glance:

- Resistant to pressure (rated for 20 bar / 290 psi continuous operation and 40 bar / 580 psi short-term)
- No cross-contamination in subsequent samples, even with demanding matrices
- Applicable for several common methods
- Quality certificate and test chromatogram for each column
- Sealed due to a clever construction with few sealing elements

Let us Refill your GPC Column

Instead of buying new columns or refilling them yourself: Simply send us your column and we will clean it, insert new filters and fill it with new material. The column comes back to you ready-to-use. Ensure a constant good quality of your columns - tested and certified (*only possible for EU states*).

Are you interested? Please contact us at service@LCTech.de or by phone +49 8082 2717-200.

Please contact us

 +49 8082 2717-0  info@LCTech.de

FREESTYLE GPC

[Automated preparation of the GPC columns via FREESTYLE GPC](#)

GPC test standard

GPC test standard for quality assurance of your GPC column.

Easy sample preparation

Fast and inexpensive access to gel permeation chromatography with
GPC uno/quattro

Brochures

[Flyer GPC columns \(pdf | 342 KB \)](#)

© 2021 LCTech GmbH

[Disclaimer](#) |

[Data Security Statement](#) |

[Terms of Service](#) |

[Sitemap](#) |