

Login 

German  

Print 

Search

*LCTech-en*

*Products*

*GPC Columns*

# GPC Columns

## Functionality Through and Through



The LCTech GPC columns are designed for easy and self-explanatory handling. Due to their robustness they can be used for every-day analysis in laboratories. The columns are available filled (with Bio-Beads) or unfilled as well as in different lengths and filling heights. They can be used with a variety of solvents.

To re-fill the columns directly on site, LCTech offers suitable kits for a simple and fast procedure. The columns can be easily taken apart, cleaned and put together again. Owing to their refined construction, the columns will be instantly tight and ready for use.





*Easy refilling of the columns*

### Further advantages at a glance:

- Resistant to pressure (rated for 20 bar / 290 psi continuous operation and 40 bar / 580 psi short-term)
- No cross-contamination in subsequent samples, even with demanding matrices
- Applicable for several common methods
- Quality certificate and test chromatogram for each column
- Sealed due to a clever construction with few sealing elements

Please contact us

+49 8082 2717-0    
info@LCTech.de

---

FREESTYLE GPC

Automated preparation of the  
GPC columns via FREESTYLE  
GPC

GPC test standard

GPC test standard for quality

assurance of your GPC column.

Easy sample preparation  
Fast and inexpensive access to gel  
permeation chromatography with  
[GPC uno/quattro](#)

Brochures

[Flyer GPC columns \( pdf | 342 KB \)](#)

© 2020 LCTech GmbH

[Disclaimer](#) |  
[Data Security Statement](#) |  
[Terms of Service](#) |  
[Sitemap](#) |