

Trade fairs cancelled. And now? LCTech presents its product news.

BE PART OF IT

Login 

German 

Print 

Search

LCTech-en

Products

Immunoaffinity columns and consumables

Aflatoxins M1

Aflatoxin M1

MYCOTOXINS

About Aflatoxin M1

Aflatoxin M1 (M = milk) is produced within the organism of animals from aflatoxin B1 after consumption with the feed and subsequently is found in milk.

Especially for this application field LCTech developed the immunoaffinity columns [AflaCLEAN M1 Select](#) and [AflaCLEAN M1 SMART](#).

The columns guarantee best results even at difficult matrices.

Both column formats are suitable for the manual processing and for fully automated processing e.g. with the LCTech robotic system [FREESTYLE SPE](#).





Sample Preparation and Analysis

Precise, simple and fast with LCTech products. In our mycotoxin brochure you find everything about mycotoxins.

Download here: [Mycotoxin brochure \(pdf|1.6 MB\)](#)

AflaCLEAN M1 Select

Immunoaffinity Columns for Clean-up of Aflatoxin M1

The AflaCLEAN M1 Select columns are available in a 3 mL polypropylene format. The quality of the columns is so high that even in difficult matrices, e.g. baby food, excellent results can be achieved.

[Read more ...](#)

AflaCLEAN M1 SMART

Immunoaffinity Columns for Clean-up of Aflatoxin M1

To speed up processing, the SMART format of the immunoaffinity columns is a good option. With a sample processing time of less than 20 minutes, the total processing time is significantly reduced in comparison with the 3 mL columns.

[Read more ...](#)

© 2021 LCTech GmbH

[Disclaimer](#) |

[Data Security Statement](#) |

[Terms of Service](#) |

[Sitemap](#) |