

Trade fairs cancelled. And now? LCTech presents its product news.

BE PART OF IT



Login 

German  

Print 

Search

LCTech-en

Products

Immunoaffinity and further Clean-up Columns

Ochratoxin A

Ochratoxin A

MYCOTOXINS

About Ochratoxin A

Ochratoxin A is a naturally occurring mycotoxin, which is highly toxic for humans and animals. Ochratoxin a is produced by various *Aspergillus* and *Penicillium* species as primary contaminant in various food and feed stuffs.

Especially for the purification of ochratoxin A in food and feed, LCTech developed the immunoaffinity columns [OtaCLEAN](#) and [OtaCLEAN SMART](#). The columns guarantee best results even at difficult matrices.

The columns possess a very high matrix tolerance and are able to bind ochratoxin A highly specific. The columns are suitable for the manual processing and for fully automated processing.



Sample Preparation and Analysis

Precise, simple and fast with LCTech products. In our mycotoxin brochure you find everything about mycotoxins.

Download here: [Mycotoxin brochure \(pdf|1.6 MB\)](#)

OtaCLEAN

Immunoaffinity Columns OtaCLEAN

The columns are available in a practical 3 mL polypropylene format. A big advantage of the columns is their long shelf life of 24 months even when stored at room temperature without compromising its high quality.

[Read more ...](#)

OtaCLEAN SMART

Immunoaffinity Columns OtaCLEAN SMART

Use the OtaCLEAN SMART columns for a shorter processing of your samples. Shelf life is 9 months cooled at 4 to 8 °C. The loading capacity is 100 ng Ochratoxin

A.

[Read more ...](#)

© 2020 LCTech GmbH

[Disclaimer](#) |

[Data Security Statement](#) |

[Terms of Service](#) |

[Sitemap](#) |