

Trade fairs cancelled. And now? LCTech presents its product news.

BE PART OF IT

Login 

German  

Print 

Search

LCTech-en

Products

Immunoaffinity columns and consumables

SPE Columns

Glass Columns

Glass Columns

Columns made of Glass for Solid Phase Extraction (SPE)

Glass columns are often needed, when plastic columns cannot be used because of interactions with the analyte or release of impurities from the synthetic material.



Available column formats

The glass columns of LCTech are available in 3 mL and 6 mL format as well as with 105 mm length and 17 mm in diameter as filled and unfilled version. The columns are filled with Florisil and sodium sulfate. Unfilled columns with 180 mm length and 17 mm in diameter are also available.

Especially for H53 method, LCTech offers the lot-certified [Elufix ready-to-use columns](#).

Automated via FREESTYLE

The LCTech glass columns (except the Version with a length of 180 mm) can be automatically processed via the [FREESTYLE SPE](#) system in combination with the inert accessories, such as stainless steel needles, glass fiber filters and caps. Subsequent evaporation or blowing off with nitrogen is possible directly within the system by

using the [EVAporation-module](#).



Adapters for 105 mm long columns, 3 mL and 6 mL columns

Please contact us

[+49 8082 2717-0](tel:+49808227170) [✉ info@LCTech.de](mailto:info@LCTech.de)

[FREESTYLE robotic system](#)

[More about the FREESTYLE robotic system and the modules in our videos.](#)

[Ready-to-use glass columns Elufix](#)

[Ready-to-use columns for H53-method](#)

[MOSH/MOAH Analytik according to DGF C-VI 22 \(20\)](#)

[Dedicated Glass Columns for MOSH/MOAH Analysis](#)

© 2021 LCTech GmbH

[Disclaimer](#) |

[Data Security Statement](#) |

[Terms of Service](#) |

[Sitemap](#) |