

Trade fairs cancelled. And now? LCTech presents its product news.

BE PART OF IT



Login

German

Print

Search

LCTech-en

Products

Immunoaffinity and further Clean-up Columns

SPE Columns

Glass Columns

Glass Columns

Columns made of Glass for Solid Phase Extraction (SPE)

Glass columns are often needed, when plastic columns cannot be used because of interactions with the analyte or release of impurities from the synthetic material.



Available column formats

The glass columns of LCTech are available in 3 mL and 6 mL format as well as with 105 mm length and 17 mm in diameter as filled and unfilled version. The columns are filled with Florisil and sodium sulfate. Unfilled columns with 180 mm length and 17 mm in diameter are also available.

Especially for H53 method, LCTech offers the lot-certified [Elufix ready-to-use columns](#).


Automated via FREESTYLE

The LCTech glass columns (except the Version with a length of 180 mm) can be automatically processed via the [FREESTYLE SPE](#) system in combination with the inert accessories, such as stainless steel needles, glass fiber filters and caps. Subsequent evaporation or blowing off with nitrogen is possible directly within the system by using the [EVAporation-module](#).



Adapters for 105 mm long columns, 3 mL and 6 mL columns

Please contact us

+49 8082 2717-0 

info@LCTech.de

FREESTYLE robotic system

[More about the FREESTYLE robotic system and the modules in our videos.](#)

Ready-to-use glass columns Elufix

[Ready-to-use columns for H53-method](#)

© 2020 LCTech GmbH

[Disclaimer](#) |

[Data Security Statement](#) |

[Terms of Service](#) |

[Sitemap](#) |