

Quality Control Certificate

Product: **Universal Column**
Product No.: 19511
Lot No.: **722388**

Storage Recommendations: Store the column at room temperature below 25°C

Description: The Universal Column is part of a 3- or 4-column setup used for the sample preparation of environmental-, food- / feed- and similar matrices with DEXTech systems from LCTech for the analysis of polychlorinated dibenzo-p-dioxins (PCDD), polychlorinated dibenzofurans (PCDF) and polychlorinated biphenyl (PCB) congeners.

Quality Control Release Inspection and Test Specification

Test Procedure: A solvent blank, spiked with quantification standard has been cleaned on a DEXTech Plus system, spiked with recovery standard, evaporated with the D-EVA and has been quantified with a HRGC/HRMS DFS from Thermo Fisher Scientific at a resolution of R > 10000.

Results Blank Value:

PCDD/F-TEQ:	0,05	pg/column
	(crit: <	0,70 pg/column)
dl-PCB-TEQ:	0,002	pg/column
	(crit: <	0,05 pg/column)
Sum Total PCB:	0,2	pg/column
	(crit: <	300 pg/column)

Results Recoveries:

PCDD/F	78	to	120	%	(crit: 70	to	120	%)
PCB	76	to	119	%	(crit: 70	to	120	%)

This is to certify that the Universal Column, Lot 722388, passed the required test specifications and is released for sale.

date: 29.07.2025

sign.: 

The company LCTech GmbH is certified according to ISO 9001

QC-Certificate - 19511 - 722388

Hazards:	<p>NOT FOR HUMAN OR DRUG USE!</p> <p>The Universal Column is designed and prepared for usage with the Alumina/Florisil Column and Carbon Column from LCTech and for laboratory use only. This product should only be used by qualified personnel familiar with its potential hazards and trained in the handling of hazardous chemicals. Due care should be exercised to prevent unnecessary human contact or ingestion, all procedures should be carried out with suitable gloves, eye protection, and clothing should be worn at all times. Waste should be disposed according to national and regional regulations.</p>
Quality Control:	<p>All ingredients are traceable to certified lots of our supplier. In addition, any ingredient with a new lot will be checked on contamination and efficiency before releasing for production. Monitoring the ongoing production, several columns are chosen at random day for analysis to check on contamination and efficiency.</p>
Quality Management:	<p>This product was produced using a Quality Management System registered to the ISO 9001:2015 (DEKRA)</p>
Documentation / Data Attached:	<p>table 1 & 2: blankvalues of PCDD/F and PCB table 3 & 4: 13C-Recoveries of PCDD/F and PCB</p>
Analytics	<p>This is to certify that the Universal Column, Lot , passed the required test specifications and is released for sale.</p>
Remarks	<p>Our suppliers maintain the highest standard of quality, however due to the high temperature necessary for several steps in the production, some small charred particles may be visible within a batch of silica or filters without any effect on the clean-up.</p>

QC-Certificate - 19511 - 722388

Results:

Lockmass check: No significant disturbances, or indicators for contaminations are detected.

Blanks: n= 6

Table 1: PCDD/F blank

	[pg/column]
2,3,7,8-TCDF	<dl
1,2,3,7,8-PeCDF	<dl
2,3,4,7,8-PeCDF	<dl
1,2,3,4,7,8-HxCDF	<dl
1,2,3,6,7,8-HxCDF	<dl
2,3,4,6,7,8-HxCDF	<dl
1,2,3,7,8,9-HxCDF	<dl
1,2,3,4,6,7,8-HpCDF	<dl
1,2,3,4,7,8,9-HpCDF	<dl
1,2,3,4,6,7,8,9-OCDF	<0,054
2,3,7,8-TCDD	<dl
1,2,3,7,8-PeCDD	<dl
1,2,3,4,7,8-HxCDD	<dl
1,2,3,6,7,8-HxCDD	<dl
1,2,3,7,8,9-HxCDD	<dl
1,2,3,4,6,7,8-HpCDD	<dl
1,2,3,4,6,7,8,9-OCDD	<dl

Table 2: PCB blank

	[pg/column]
PCB-#28	0,19
PCB-#52	<0,144
PCB-#101	<dl
PCB-#153	<dl
PCB-#138	<dl
PCB-#180	<dl
PCB-#81	<dl
PCB-#77	<0,045
PCB-#126	0,0163
PCB-#169	<0,027
PCB-#123	0,02
PCB-#118	<dl
PCB-#114	0,02
PCB-#105	<dl
PCB-#167	<0,027
PCB-#156	<dl
PCB-#157	0,02
PCB-#189	0,015

PCDD/F TEQ (2005)	[pg/column]
lower bound	0
upper bound	0,05

PCB-TEQ	[pg/column]
lower bound	0,002
upper bound	0,002
Sum DIN	0,2

QC-Certificate - 19511 - 722388

Table 3: PCDD/F recoveries

		[%]	RSD [%]
PCDD/F 13C Recoveries [%]	2,3,7,8-TCDF	98	3
	1,2,3,7,8-PeCDF	92	2
	2,3,4,7,8-PeCDF	90	4
	1,2,3,4,7,8-HxCDF	90	10
	1,2,3,6,7,8-HxCDF	91	10
	2,3,4,6,7,8-HxCDF	86	10
	1,2,3,7,8,9-HxCDF	87	10
	1,2,3,4,6,7,8-HpCDF	86	7
	1,2,3,4,7,8,9-HpCDF	120	10
	1,2,3,4,6,7,8,9-OCDF	111	0
	2,3,7,8-TCDD	96	9
	1,2,3,7,8-PeCDD	89	2
	1,2,3,4,7,8-HxCDD	88	10
	1,2,3,6,7,8-HxCDD	78	8
	1,2,3,7,8,9-HxCDD	92	8
	1,2,3,4,6,7,8-HpCDD	120	9
	1,2,3,4,6,7,8,9-OCDD	116	12

Table 4: PCB recoveries

		[%]	RSD [%]
PCB 13C Recoveries [%]	PCB-#28	97	5
	PCB-#52	90	9
	PCB-#101	97	4
	PCB-#153	81	9
	PCB-#138	83	8
	PCB-#180	89	5
	PCB-#81	116	3
	PCB-#77	111	8
	PCB-#126	119	2
	PCB-#169	112	7
	PCB-#123	82	12
	PCB-#118	76	16
	PCB-#114	97	8
	PCB-#105	81	20
	PCB-#167	79	10
	PCB-#156	85	6
	PCB-#157	77	9
	PCB-#189	76	10