

Quality Control Certificate

Product: **Smart Column**
 Product No.: 19513
 Lot No.: **719528**

Storage Recommendations: Store the column at room temperature below 25°C

Description: The Smart Column is part of a 3-column setup used for the sample preparation of environmental-, food- / feed- and similar matrices with DEXTech systems from LCTech for the analysis of polychlorinated dibenzo-p-dioxins (PCDD), polychlorinated dibenzofurans (PCDF) and polychlorinated biphenyl (PCB) congeners.

Quality Control Release Inspection and Test Specification

Test Procedure: A solvent blank, spiked with quantification standard has been cleaned on a DEXTech Plus system, spiked with recovery standard, evaporated with the D-EVA and has been quantified with a HRGC/HRMS DFS from Thermo Fisher Scientific at a resolution of R > 10000.

Results Blank Value:	PCDD/F-TEQ:	0,12	pg/column
		(crit: <	0,7 pg/column)
	dl-PCB-TEQ:	0,0257	pg/column
		(crit: <	0,05 pg/column)
	Sum Total PCB:	15	pg/column
		(crit: <	300 pg/column)

Results Recoveries:	PCDD/F	77	to	109	%	(crit: 70	to	120	%)
	PCB	70	to	98	%	(crit: 70	to	120	%)

This is to certify that the Smart Column, Lot 719528, passed the required test specifications and is released for sale.

date: 25.01.2024 sign.: *M. Brachis*

The company LCTech GmbH is certified according to ISO 9001



QC-Certificate - 19513 - 719528

Hazards:	<p>NOT FOR HUMAN OR DRUG USE!</p> <p>The Smart Column is designed and prepared for usage with the Alumina/Florisil Column and Carbon Column from LCTech and for laboratory use only. This product should only be used by qualified personnel familiar with its potential hazards and trained in the handling of hazardous chemicals. Due care should be exercised to prevent unnecessary human contact or ingestion, all procedures should be carried out with suitable gloves, eye protection, and clothing should be worn at all times. Waste should be disposed according to national and regional regulations.</p>
Quality Control:	<p>All ingredients are traceable to certified lots of our supplier. In addition, any ingredient with a new lot will be checked on contamination and efficiency before releasing for production. Monitoring the ongoing production, several columns are chosen at random day for analysis to check on contamination and efficiency.</p>
Quality Management:	<p>This product was produced using a Quality Management System registered to the ISO 9001:2015 (DEKRA)</p>
Documentation / Data Attached:	<p>table 1 & 2: blankvalues of PCDD/F and PCB table 3 & 4: 13C-Recoveries of PCDD/F and PCB</p>
Analytics	<p>This is to certify that the Smart Column, Lot , passed the required test specifications and is released for sale.</p>
Remarks	<p>Our suppliers maintain the highest standard of quality, however due to the high temperature necessary for several steps in the production, some small charred particles may be visible within a batch of silica or filters without any effect on the clean-up.</p>



QC-Certificate - 19513 - 719528

Results:

Lockmass check: No significant disturbances, or indicators for contaminations are detected.

Blanks: n= 6

Table 1: PCDD/F blank

	[pg/column]
2,3,7,8-TCDF	<dl
1,2,3,7,8-PeCDF	0,09
2,3,4,7,8-PeCDF	0,12
1,2,3,4,7,8-HxCDF	0,076
1,2,3,6,7,8-HxCDF	0,033
2,3,4,6,7,8-HxCDF	<0,045
1,2,3,7,8,9-HxCDF	0,09
1,2,3,4,6,7,8-HpCDF	<0,063
1,2,3,4,7,8,9-HpCDF	0,018
1,2,3,4,6,7,8,9-OCDF	0,06
2,3,7,8-TCDD	<0,036
1,2,3,7,8-PeCDD	<dl
1,2,3,4,7,8-HxCDD	<0,027
1,2,3,6,7,8-HxCDD	<0,108
1,2,3,7,8,9-HxCDD	0,068
1,2,3,4,6,7,8-HpCDD	<0,09
1,2,3,4,6,7,8,9-OCDD	0,46

Table 2: PCB blank

	[pg/column]
PCB-#28	3,13
PCB-#52	4,06
PCB-#101	2,97
PCB-#153	2,09
PCB-#138	2,04
PCB-#180	0,695
PCB-#81	0,06
PCB-#77	0,283
PCB-#126	0,2418
PCB-#169	0,045
PCB-#123	0,36
PCB-#118	1,35
PCB-#114	0,252
PCB-#105	0,56
PCB-#167	0,194
PCB-#156	0,248
PCB-#157	0,15
PCB-#189	0,345

PCDD/F TEQ (2005)	[pg/column]
lower bound	0,1
upper bound	0,12

PCB-TEQ	[pg/column]
lower bound	0,0257
upper bound	0,0257
Sum DIN	15

QC-Certificate - 19513 - 719528

Table 3: PCDD/F recoveries

	[%]	RSD [%]	
PCDD/F 13C Recoveries [%]	2,3,7,8-TCDF	88	5
	1,2,3,7,8-PeCDF	77	10
	2,3,4,7,8-PeCDF	80	12
	1,2,3,4,7,8-HxCDF	100	7
	1,2,3,6,7,8-HxCDF	109	6
	2,3,4,6,7,8-HxCDF	103	3
	1,2,3,7,8,9-HxCDF	96	7
	1,2,3,4,6,7,8-HpCDF	99	5
	1,2,3,4,7,8,9-HpCDF	90	8
	1,2,3,4,6,7,8,9-OCDF	89	9
	2,3,7,8-TCDD	81	8
	1,2,3,7,8-PeCDD	80	14
	1,2,3,4,7,8-HxCDD	105	6
	1,2,3,6,7,8-HxCDD	87	4
	1,2,3,7,8,9-HxCDD	105	4
	1,2,3,4,6,7,8-HpCDD	93	3
	1,2,3,4,6,7,8,9-OCDD	89	3

Table 4: PCB recoveries

	[%]	RSD [%]	
PCB 13C Recoveries [%]	PCB-#28	86	10
	PCB-#52	79	8
	PCB-#101	93	7
	PCB-#153	90	6
	PCB-#138	87	5
	PCB-#180	98	9
	PCB-#81	77	8
	PCB-#77	79	4
	PCB-#126	70	11
	PCB-#169	71	19
	PCB-#123	87	6
	PCB-#118	83	8
	PCB-#114	88	6
	PCB-#105	86	6
	PCB-#167	84	11
	PCB-#156	96	12
	PCB-#157	96	11
	PCB-#189	90	9