

Quality Control Certificate

Product: **EVOLUTION Universal Column**
Product No.: 20085
Lot No.: **718969**

Storage Recommendations: Store the column at room temperature below 25°C

Description: The EVOLUTION Universal Column is part of a 3- or 4-column setup used for the sample preparation of environmental-, food- / feed- and similar matrices with DEXTech systems from LCTech for the analysis of polychlorinated dibenzo-p-dioxins (PCDD), polychlorinated dibenzofurans (PCDF) and polychlorinated biphenyl (PCB) congeners.

Quality Control Release Inspection and Test Specification

Test Procedure: A solvent blank, spiked with quantification standard has been cleaned on a DEXTech Plus system, spiked with recovery standard, evaporated with the D-EVA and has been quantified with a HRGC/HRMS DFS from Thermo Fisher Scientific at a resolution of R > 10000.

Results Blank Value:

PCDD/F-TEQ:	0,15	pg/column
	(crit: <	0,7 pg/column)
dl-PCB-TEQ:	0,0063	pg/column
	(crit: <	0,05 pg/column)
Sum Total PCB:	16,9	pg/column
	(crit: <	300 pg/column)

Results Recoveries:

PCDD/F	83	to	117	%	(crit: 70	to	120	%)
PCB	72	to	100	%	(crit: 70	to	120	%)

This is to certify that the EVOLUTION Universal Column, Lot 718969, passed the required test specifications and is released for sale.

date: 23.08.2023 sign.: T. Kehmeier

The company LCTech GmbH is certified according to ISO 9001



QC-Certificate - 20085 - 718969

Hazards:	<p>NOT FOR HUMAN OR DRUG USE!</p> <p>The EVOLUTION Universal Column is designed and prepared for usage with the Alumina/Florisil Column and Carbon Column from LCTech and for laboratory use only. This product should only be used by qualified personnel familiar with its potential hazards and trained in the handling of hazardous chemicals. Due care should be exercised to prevent unnecessary human contact or ingestion, all procedures should be carried out with suitable gloves, eye protection, and clothing should be worn at all times. Waste should be disposed according to national and regional regulations.</p>
Quality Control:	<p>All ingredients are traceable to certified lots of our supplier. In addition, any ingredient with a new lot will be checked on contamination and efficiency before releasing for production. Monitoring the ongoing production, several columns are chosen at random day for analysis to check on contamination and efficiency.</p>
Quality Management:	<p>This product was produced using a Quality Management System registered to the ISO 9001:2015 (DEKRA)</p>
Documentation / Data Attached:	<p>table 1 & 2: blankvalues of PCDD/F and PCB table 3 & 4: 13C-Recoveries of PCDD/F and PCB</p>
Analytics	<p>This is to certify that the EVOLUTION Universal Column, Lot , passed the required test specifications and is released for sale.</p>
Remarks	<p>Our suppliers maintain the highest standard of quality, however due to the high temperature necessary for several steps in the production, some small charred particles may be visible within a batch of silica or filters without any effect on the clean-up.</p>



QC-Certificate - 20085 - 718969

Results:

Lockmass check: No significant disturbances, or indicators for contaminations are detected.

Blanks: n= 6

Table 1: PCDD/F blank

	[pg/column]
2,3,7,8-TCDF	0,04
1,2,3,7,8-PeCDF	0,09
2,3,4,7,8-PeCDF	<0,081
1,2,3,4,7,8-HxCDF	0,037
1,2,3,6,7,8-HxCDF	<0,018
2,3,4,6,7,8-HxCDF	<dl
1,2,3,7,8,9-HxCDF	<0,045
1,2,3,4,6,7,8-HpCDF	0,11
1,2,3,4,7,8,9-HpCDF	0,018
1,2,3,4,6,7,8,9-OCDF	0,06
2,3,7,8-TCDD	<0,036
1,2,3,7,8-PeCDD	0,09
1,2,3,4,7,8-HxCDD	<0,027
1,2,3,6,7,8-HxCDD	<0,108
1,2,3,7,8,9-HxCDD	<0,027
1,2,3,4,6,7,8-HpCDD	0,12
1,2,3,4,6,7,8,9-OCDD	1,15

Table 2: PCB blank

	[pg/column]
PCB-#28	6,35
PCB-#52	5,35
PCB-#101	1,79
PCB-#153	1,26
PCB-#138	1,19
PCB-#180	0,927
PCB-#81	0,11
PCB-#77	0,168
PCB-#126	0,0483
PCB-#169	0,041
PCB-#123	0,53
PCB-#118	1,29
PCB-#114	0,618
PCB-#105	0,74
PCB-#167	0,543
PCB-#156	0,653
PCB-#157	0,55
PCB-#189	0,738

PCDD/F TEQ (2005)	[pg/column]
lower bound	0,15
upper bound	0,15

PCB-TEQ	[pg/column]
lower bound	0,0063
upper bound	0,0063
Sum DIN	16,9

Table 3: PCDD/F recoveries

	[%]	RSD [%]	
PCDD/F 13C Recoveries [%]	2,3,7,8-TCDF	95	7
	1,2,3,7,8-PeCDF	93	6
	2,3,4,7,8-PeCDF	95	3
	1,2,3,4,7,8-HxCDF	108	7
	1,2,3,6,7,8-HxCDF	117	8
	2,3,4,6,7,8-HxCDF	107	7
	1,2,3,7,8,9-HxCDF	111	6
	1,2,3,4,6,7,8-HpCDF	105	6
	1,2,3,4,7,8,9-HpCDF	98	3
	1,2,3,4,6,7,8,9-OCDF	89	3
	2,3,7,8-TCDD	91	4
	1,2,3,7,8-PeCDD	91	4
	1,2,3,4,7,8-HxCDD	106	5
	1,2,3,6,7,8-HxCDD	94	4
	1,2,3,7,8,9-HxCDD	110	4
	1,2,3,4,6,7,8-HpCDD	98	3
	1,2,3,4,6,7,8,9-OCDD	83	5

Table 4: PCB recoveries

	[%]	RSD [%]	
PCB 13C Recoveries [%]	PCB-#28	90	1
	PCB-#52	100	6
	PCB-#101	95	5
	PCB-#153	96	3
	PCB-#138	91	3
	PCB-#180	84	4
	PCB-#81	82	4
	PCB-#77	81	5
	PCB-#126	72	4
	PCB-#169	72	5
	PCB-#123	79	3
	PCB-#118	74	2
	PCB-#114	81	3
	PCB-#105	75	4
	PCB-#167	86	2
	PCB-#156	74	2
	PCB-#157	73	4
	PCB-#189	74	7