Approved and Ready-to-Use

"Just click it": Outstanding results with proven chemistry for analysis of Dioxine/PCB and other POPs. Take the columns out of the package and click them into the column tower with a simple one-handed operation. No screws, no tools, or additional steps needed. Depending on your sample and chosen method, glass or plastic columns can be inserted into the 3-columns set-up with equivalent results.



Click into any DEXTech System

Carefully single wrapped in polybag packing



Top view of 100 plastic columns



Dioxin Glass Columns

- Proven glass body
- Excellent results
- Very high loading capacity



EVOLUTION Plastic Columns

- Phthalate free (Chemical inert PP)
- Unbreakable
- Very high loading capacity
- Absolutely equivalent results to glass version

Column Overview

glass columns

	Name	P/N	Material	pcs	Food	Feed	Environmental	Analytes
Position 1	Universal	19511	Glass	25/box	up to 5 g fat	х	x	PCDD/F, PCBs, PBDEs, PCNs
	Universal EVOLUTION	20085	Plastic	100/box	up to 5 g fat	х	х	PCDD/F, PCBs, PCN
	Standard	19512	Glass	25/box	up to 5 g fat	х		PCDD/F, PCBs, PBDEs, PCNs
	PCB209	19871	Glass	25/box	х		x	all 209 PCBs
	PCB209 EVOLUTION	20325	Plastic	100/box	x		х	all 209 PCBs
	SMART	19513	Glass	25/box	up to 1.5 g fat	Х	х	PCDD/F, PCBs, PBDEs, PCNs, CPs
Pos. 2	Alox	15433	Glass	25/box	up to 5 g fat	х	x	PCDD/F, PCBs, PBDEs, PCNs
	Alox EVOLUTION	20087	Plastic	100/box	up to 5 g fat	х	x	PCDD/F, PCBs, PBDEs, PCNs
Pos. 3	Carbon column	20777	stainless steel	25/box	up to 5 g fat	Х	х	PCDD/F, PCBs, PBDEs, PCNs